

Product datasheet for **SC401045**

Mir449b Mouse MicroRNA Expression Plasmid (MI0005547)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: AGACUCGGGUAGGCAGUGUUGUAGCUGGCUGCGUUGGUCAGGCCAGCAGCCACAGCUACCCUGCCACUCCUUCUGGC
Mature Sequence: AGGCAGUGUUGUAGCUGGC

Full insert: >MI0005547 AATTTTGTCTTATAGTTTAGGAGACATTAAACAAGGTTTCTCTATAAATCCTCTTCAGTTACTAAGAT
 TAGAGTTCTGGGTTGGCAAGAAGGTTCAAGTAAAGGCACTTACTGCCATGCTATCCCAGGAATTCACATG
 GTGGAAGGAGAAAAGTGGCTCCCTCAGAGTTCCCTTGCTCTCCGTATCCACGCCACCACACTCCAGCAGC
 TCTGCAGCAGTGTCTTCTGTTAGAACTTAGCCTCAGCCACTCCACAGACTGTGAAAATCCGAACTTGC
 CAGACTCGGGTAGGCAGTGTGTTAGCTGGCTGCGTTGGGTCAGGCCAGCAGCCACAGCTACCCTGCCAC
 TTCCTTCTGGCTGGGTTCTTCTGTTAATGAGGTGCGCTGCTTGCTGTGTGTGATGGCTTGGCAGTGT
 ATTGTTAGCTGGTTGAGTATGTGAGCGGCACCAGCTAACATGCGACTGCTCTCCTATTGCACACACAATC
 ACCACCAGCATGCATTTACAGCTCATGTATGTGCTTTGCCATGAAGAGGGGTACTATTATGATAACAGCT
 ATTAGTGGCTTGGGGTGACACAAGGACATATTTATAGCAGTGAAAGTACTAACCATTGTTACTTTAAAAA
 CTCAAGTAATATTTTCTGAGTCCTAGCTTG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: [MI0001649](#), [MI0005547](#)
 Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0005547](#)


[View online »](#)