

Product datasheet for SC401047

Mir450a-1 Mouse MicroRNA Expression Plasmid (MI0001653)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GAGAGAUACUGAGCUGUUUUUGCGAUGUGUCCUAAUAUGUGCUAUAUUUAUUAUUGGGAACAUUUUGCAUAAAUAGCUUUGUGUCAAUACA
Mature Sequence:	UUUUGCGAUGUGUCCUAAUAU
Full insert:	>MI0001653 TTCTCCGTTTTAATTTTACATACATGGGATCAATATACAGATATAAGGCTAAGAGAACAGATCAAAATAT TATAGTTATCCTGAGTTGTGCTTTAATGTTATTGAACAAATATGCTGGGACTGATGCTGTGGTTGAATG ATAGAGGAAGGAGAAAGGCAAACCAAGATGCTGAACTATTTTTCGATGTGTTCTAATATGTACTATGA GTATATTGGGGATGCTTTGCATTCATGGTTCTGCATCCATCATATGGACAGTCACCTCCAAGCTGTAGCA CCACCACTGACCAAAGAGAGATACTGAGCTGTTTTTTCGATGTGTTCTAATATGTGCTATAATTATATT GGGAACATTTTGCATAAATAGCTTTGTGTCAATACACAGGCAAAAGGTTGGAAGAAAGAAGCTGCAAGTC TAACCATAGTTGTTTATGGTTTGTATTAATAGGAATCAGATACAGAATTATTTTGCAGTATGTTCC TGAATACATAACATAAGCGTATTGGGAACATTTTGCATGCATAATTTTGTATCAATATATGGGGGAAGA CTAGGAAGAAACATGACAAAACCTGAACAGAACTTATGTCTCTGGTATTACAGGTGATGTTTATTGTTGT CTTTATATTTTATTGTTTTTTTTTACCAATTATTGTGTAACAAAATGTTT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	MI0001653, MI0003537, MI0004705 Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0001653


[View online »](#)