

Product datasheet for **SC401066**

Mir466d Mouse MicroRNA Expression Plasmid (MI0005546)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: CAUGUGUGUUUGUGUGUGCGUACAUGUACAUGUGUGUUAUGAAUUAUCAUAUACAUAACACGCACACAUAGAUACGCACGCACACACACACAGG
Mature Sequence: UGUGUGUGCGUACAUGUACAUG

Full insert: >MI0005546 TGCTTGTGTATATGTATGTGTGCATGAACATGTGTATATGAATATACATATACATACACACATACCA
TACAAACATGCATACAAACACACAGAAAATGGATACTCATTGTTCTATTGAATGATATGTGTCTTAGAGT
CTTAGAGATTCTGTATGACTGGGCCAGGGAGAGGTCACACTCACCATGTGACGCCCTTGCCAGAGTCAG
GGTGATCAGGAGTGCAAATTCATGTTAACTGCGTGATGAGGTAGAGATCAGCCTGTCATTCATGTGATGG
CCCCTGTGTATGTGCATGTGTGTTGTGTGTGCGTACATGTACATGTGTGTATATGAATATACATATACA
TACACGCACACATAGATACGCACGCACACACACACAGGGATGGCACTCATCAATCCATTTGAATGATA
AGTGACTTAGAGTGAAAATGCTTACTAGAGTTAAGAGTATCAGGAGTTCAAATTCATACTTATCCTCATG
AGAATTTAGAGATCAGCCTGTGATTTATGTGATTCTCCATGTATGTGCATGTGTGTATAGTTGTGTGTG
CATGTTTCATGTCTATATTTGAATATACATAACATACACACACACAGTATAAACGCAAGCACACATACACAG
ATACAGGAATGGCATTGCGCAATCCCTTTGAATGATAAGTC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: [MI0004523](#)
Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0005546](#)



[View online »](#)