

Product datasheet for SC401145

Mir669h Mouse MicroRNA Expression Plasmid (MI0006289)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: UGUGAUGCCCAAAAAUAUGCAUGGGUGUAUAGUUGAGUGCAUUAUUAUUCUCUAUAUGAAUAUGCAUUAACACACAUGCACACACACACACUCACACAC

Mature Sequence: AUGCAUGGGUGUAUAGUUGAGUGC

Full insert: >MI0006289 GAGTTAGAGATCAGCCTGGAATTCATTTGATTGTCCTTGTGTATGTGTATATATGTGTAGGTAAGTGTGC
 TTGTACATGCATGTATATGAACATACACACACATACATATACACATAAGAATGGCAGTCATAGATCTG
 TTTGATTTATATGTGCTTAGAGTTTATAGAGTGTATAATTATGGGTCAGTGTTACATCTAACTCACCAG
 GTGAAAACCACTTGCCAGAGTTAGAAGGATCAGGAATTTAAATTCATGCTTAACTTCCTGATAAGTTTGT
 GATCATCCTGGCATTATGTGATGCCCAAAAAATATGCATGGGTGTATAGTTGAGTGCATATTCATTT
 CTCTATATGAATATGCATATACACACATGCACACACACACACTCACACACAGGAATGACACTCGTTTA
 TCCTTTTGAATGGTATGTGACTCAGAGCTTTAGAGCTTATGCAGTACTGGGCCAGAGCTATGTCTCAGTT
 AGGAGGTGAATACTCTTGCCAGCGTCAGGAGGTCCTGGAGTTCAAATTCATGTTTCACTTCATGATGAGT
 TAGAGATCAGCATGCGATTATGTGATGCCCTCCCCATGTATGTGCATGGGTTCACTTATGTGTGCA
 TGTAATGTATGTATGAATGTTACATACATATACATAAACACACACACTGGCATTACCTACACGCA
 CATTAGTGGCTCTG

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect
 the mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0006289](#)



[View online »](#)