

## Product datasheet for SC401204

### Mir715 Mouse MicroRNA Expression Plasmid (MI0004700)

#### Product data:

**Product** MicroRNA Expression Plasmids

**Type:**

**Species:** Mouse

**Vector:** pCMV-MIR (PCMVMIR)

**Precursor:** GGCUGGGGAGAGGGCUCCGUGCACACCCCGCGUGCGCGGUACUUUCCUCCCUCCUGAGGGCCGCCGUGCGGGACGGGGUGUGGGUAGGCAACGGUGGGCUCCG

**Mature** CUCCGUGCACACCCCGCGUG

**Sequence:**

**Full insert:** >MI0004700 TCGTCGGGGTTTTGGGTCCGTCGCCCTCAGTGAGAAAGTTTCCTTCTCTAGCTATCTTCCGAAAGGG  
 TGC GGCTTCTTACGGTCTTGAGGGGTCTCTCCGAATGGTCCCCTCGAGGGCTCGCCCCCTGACCGCT  
 CCCGTGCGCGCAGCGTTTGTCTCTCGTCTACACGGCCAGCGCTCCCGCTCCGAGTTCGGGGAGGGA  
 TCACGCGGGCAGAGCTGTCTGTCTCTGCGTCCGCTGCGGAGCCATGTGGCTCCGTTGTGTGGTGG  
 TGGCTGGGGAGAGGGTCCGTGCACACCCCGCGTGC GCGGTACTTTCTCCCTCCTGAGGGCCGCCGT  
 GCGGGACGGGTGTGGGTAGGCAACGGTGGGTCCCGGTCCCAACCGTCTTCCCGTGCCTACCCGTG  
 CCTTCGTCGCGTGCCTCTCGCTCGCGTCCACGACTTTGGCCGCTCCCGCGACGCGGGCCTGCGCC  
 GCGGTGGTGGGTGTGTGTCTCGGGTGTGTGGTTGTGTGCGCTCGCCCCCTTCCCGCGGCAG  
 CGTTCACGCTGGCGAAGTCGCGGGAGTCTCTCCCTCTTCCGGGTGAGAGGGTCCGTGTCTGG  
 CGTTGATTGATCTCGCTCTCGGGACGGGACCGTTCTGTGGGAGAACGGCTGTTGG

**OTI** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The  
**Disclaimer:** flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms.  
 This should not affect the function of the mature miRNA.

**E. coli** Kanamycin (25 ug/mL)

**Selection:**

**Mammalian** Neomycin

**Cell**

**Selection:**

**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

**Also** Both mature and minor (star) miRNA can be over expressed from the construct.

**Contains:**

**RefSeq:** [MI0004700](#)



[View online »](#)