

Product datasheet for SC401225

Mir7b Mouse MicroRNA Expression Plasmid (MI0000730)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	AGGAGCGGAGUACGUGAGCCAGUGCUAUGUGGAAGACUUGUGAUUUUGUUGUUCUGAUUAUGAUUAUGACAACAAGUCACAGCCAGCCUCAUAGCGUGGACUCCUAU
Mature Sequence:	UGGAAGACUUGUGAUUUUGUUGU
Full insert:	>MI0000730 CTTGAAGTTCTCATCCTCCTGCCTCTACTCCCTGCGTGCTGGGGTGACAAGTGTGTGCTCCTCGTGCCTG GAACAATAAATCCCACGTTGGTGATGTGCCAGGGTGTGAGTTGGGGGCAGGTAGAGAATGAATCTTGCC TGTGTCTCAGGGGATGGAATGGAAAAGGGATGTCTAGGGAGGGTAAGTCAGTGAGTGGGGTTTGTATGC TCAATCCTCTATGGAGAGAGAGAGAGAAGCACTTGAGGGTGTGAACCTGGCGGAGGAGGATGAGGACGA GGAGCGGAGTACGTGAGCCAGTGCTATGTGGAAGACTTGTGATTTTGTGTTCTGATATGATATGACAAC AAGTCACAGCCAGCCTCATAGCGTGGACTCCTATCACCTTTCTCTGAAATACACAAACAGGTAGGGGCA GGAAACTGGGACAGGAGGGGGTTAAGGGAAGTCTAGCACGGACAGGAGGGTGTCTAGCAAGGTCTCGG GCCCTGGGAGAGTGCAGGAGTGTTCCTATGCTATTGGCCTGTCACTCAAACACATCTACATGATAGGAT CAGAGACACAACCTCAGGGCCCTGGGTACAGCACCGTGCTTCTCTGGAGACTCAAATAGGGAATGGTG AGCCAGTGGAGTAACTCAGAGGGAAGAGTGCTTTGAA
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The
Disclaimer:	flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000730
Synonyms:	Mirn7b mmu-mir-7b
Locus ID:	723883


[View online »](#)