

## Product datasheet for **SC401248**

### Mir93 Mouse MicroRNA Expression Plasmid (MI0000581)

#### Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	AGUCAUGGGGGCUCCAAGUGCUGUUCGUGCAGGUAGUGUAAUUACCUGACCUACUGCUGAGCUAGCACUUCCCGAGCCCCAGGACA
Mature Sequence:	CAAAGUGCUGUUCGUGCAGGUAG
Full insert:	>MI0000581 GTGAGTGGTGGGTCCTGTCAGCTGGAAAGCTGACCCCTGCCTGCCACCTCACCTAATGACCCTCAAGCC GCCTCCTTCCCTCCTACCAGCCCTGCTGGGACTAAAGTGTGACAGTGCAGATAGTGGTCTCTGTGC TACCGCACTGTGGTACTTGTCTGCTCCAGCAGGGCACGTGCAACGCCCATGGAGGAAAGGCTGCTTGTCT GCTTGAATCCATGAGGGCTAGGACTTGAATTCCTGGTGTCTAGATATGATCTTAGAAATCTCCATTCTG TCCTCCCCTCTTGGACCTTAGTCATGGGGCTCCAAAGTGTGTTTCGTGCAGGTAGTGAATTACCTGAC CTACTGCTGAGCTAGCACTTCCCAGCCCCAGGACACAACCTCTGATAGTGGCCCTAGCTGAGCGCCCC AGGCTCCATTTGTTGATGGGGAGCCTAGTGGAGTTCAGAAGGGTCTGGTCTCCCTCACAGGACAGCTGA ACACCCACCGGGACTGGCCAGTGTGAGAGCGGAGACTTGGGCAATTGCTGGACGCTGCCCTGGGCATT GCACTTGTCTCGGTCTGACAGTGCCGGCCCAACACTGCAGATGTGGGGGGCAGGGAGGAGTCCCCATTT TTGTTCAATGAGTAGGCAAATGTGGCTTGCCAGGCTTGAAGCTGGAGCAAAGA
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	<a href="#">MI0000407</a> , <a href="#">MI0000581</a> , <a href="#">MI0000689</a> Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	<a href="#">MI0000581</a>
Synonyms:	Mir93 mmu-mir-93
Locus ID:	723885



[View online »](#)