

Product datasheet for SC401255

MIR1184-3 Human MicroRNA Expression Plasmid (MI0015972)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CUUGCAGAACGAGGUGAAGGAGGUGGUUCUGCUCAGCAGUCAACAGUGGCCACAUCUCCACCUGCAGCGACUUGAUGGCUUCCGUGUCCUUUUCGUGGG
Mature Sequence:	CCUGCAGCGACUUGAUGGCUUCC
Full insert:	>MI0015972 CACCCCAGCCCGGGCCCGGGGCGACGCCCGCCCTACCGGCCGCGTGCTTCTCCGAACTCCGGCTCGGC GGCGCCCTCTCCCGCCGCCCTGGGCGCCTTCTCGGACGTGCTGGTCCGCTGCGAGGTGTCCCGCGTGCTG CTGCTGCTCCTCTGCAACCACCGCCGCCAAGCTGCTGCCGAGCAGCCCAGACCCTGGAGAAGTACT CCTGGGAGGCTTTTGACAGCCACGGGCAGGAGAGCAGCGGCCAGCTTCCCGAGGAGCTTTTCTGCTGCT CCAGTCTTTGGTCATGGCTACCCACGAAAAGGACACGGAAGCCATCAAGTCGCTGCAGGTGGAGATGTGG CCACTGTTGACTGCTGAGCAGAACCACCTCCTTACCTCGTTCTGCAAGAAACCATCTCCCCCTCAGGAC AGGGAGTCTGATCCATCCCATTACCCAGTGACTTCTTTTTGCCAGGCTGGACTTTTTGCATCAGTCA CGTTAACAGATGACTTTGCCTGTTACCAAACCTCATGCATCCACGTTTGCCTGTTGGGAGGAATAAAAA GACATCGTTCCCGCTTCTGCGTTTTGTTATTCTACTGCCGCCATAGGAATTATTCGTTGGCTGAACGT TACCAGCACCCCGAGAACACATTTTGATAGAATCAGAGTAGAGGACATGGCTGTCTTCTAAAAAGCCAC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0015972
Synonyms:	MIR1184 MIRN1184 hsa-mir-1184 hsa-mir-1184-1
Locus ID:	100302111



[View online »](#)