

Product datasheet for SC401294

MIR2682 Human MicroRNA Expression Plasmid (MI0012063)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: ACCUUCUGAAAGAGGUUGGGGAGGCAGUGACUGUUCAGACGUCCAUCUCUUUGGGACGCCUCUUCAGCGUGUCUCCUGCCUCUGCCUUUAGGACGAGUC

Mature Sequence: CAGGCAGUGACUGUUCAGACGUC

Full insert: >MI0012063 GAAAGGGCGTGTCCGTGGGCGAAGGTGCCCCGCAAGCCTGCGGAAGCGCGTGACTTCGCGAATCCCGAG
CCCCGGGCGGTGCCGGAGCCCCACCCCGGCCGCGCTGCTCCGCAACCAGACTTGGCGGCCGAGAT
GGGGAGGGGGCGGGGGCCGGGGGCTGGGTCCCCGCGGAAGTCCAGATCGGAAGTGTGGCGGGTCTCTC
AGGAACGCATATTTTCTTTTCGGATTCTCCGAGAAGGGGTGGAGAAACAGGATGGGAACCGGATTGGTTG
CAGTGTGGCACAGACCTGGAACCTTCCTGAAAGAGGTTGGGGCAGGCAGTGACTGTTTCAGACGTCCAATC
TCTTTGGGACGCTCTTCAGCGCTGTCTTCCCTGCCTCTGCCTTTAGGACGAGTCTCAAACACCAACAAA
CGCAAGGCACGTCCCCCTCTCAGGTAGCTGCGAAGGCGTGGGAGCAAGAGAGAGCCTCTCCCTCCTTC
TACTTAGAGTTTTTGTGTTGTTGTTGTTGTTGTTGTTTCAAATATTAATGGACGATAAAACAACACGGAT
TTGTCACCTCTAGCGAAGAGGAAGCCACCGCCACCACCCCTCAAGCATGCCTGGGCACTCTCTGGGGA
CAATCGGACGTGGCAGGAGTTGGCATGCCAAGATGAAGCTGGGAAGGATGACCCAGAGTGGAAGGCGTG
GGCTGGGGCT

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms.
This should not affect the function of the mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0012063](#)

Locus ID: 100616452



[View online »](#)