

Product datasheet for SC401295

MIR2861 Human MicroRNA Expression Plasmid (MI0013006)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GGCGCCUCUGCAGCUCCGGCUCUCCCCUGGCCUCUCGGGAACUACAAGUCCAGGGGGCCUGGCCGUGGGCGGCGGCGGAAGAGGCGGGG
Mature Sequence:	GGGGCCUGGCCGUGGGCGG
Full insert:	>MI0013006 CGGGCCCGGGTCCGGCCAGGGAGCTCCAGGCCGCGGCTCGGGGGCTGACCCCCCGCCGCCCCCGCG CAGCCCGGCGGAGGCGGCGCCGTCGGAGGAGGAGGAGCAGGGCACTGCGGCTGGCCCCGGATCGGGCGCC AGAGCGGCAGCAGCTTCGGCAGCAGCGGCGGCCCGGGCCCATGCAGCGGGACGCGCCACCCCGAGCCCCA GCTCCGGCGCCCCGGCTCCCCGCGCCCCGATCGGGGCGCGCGCTAGTAGTGGCGGCGGCGGAGGCGGGG GCAGCGGCGGCGGCGGCGGAGGCGCCTCTGCAGCTCCGGCTCCCCCTGGCCTCTCGGGAACATAAGTCC CAGGGGGCCTGGCGGTGGGCGGCGGCGGAAGAGGCGGGTTCGGCGCCGCGAGGCCGGAAGTGGCCGTGG AGGCGGAAGTGGCGGCGCGCGGAGGGGCTGGAGTGCGGCGGCGGCGGACCCGGAGCAGGAGCGGCGG CAGCAGCGACTGGGGGCGGCGGCGCGCTTGAGGCGGCCATGGCAAAGCAGTACGACTCGGTGGAGTG CCCTTTTGTGATGAAGTTTCAAATACGAGAAGCTCGCCAAGATCGGCCAAGGCACCTTCGGTAAGGC TGGGCCCTCGGGCGGCGGAGCCCTGGCCCTGCACCCCTAGGGCCGACGTCGGGATGCCC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0013006
Synonyms:	BMND15 MIRN2861
Locus ID:	100422910


[View online »](#)