

Product datasheet for **SC401356**

MIR3158-1 Human MicroRNA Expression Plasmid (MI0014186)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: AUUCAGGCCGGUCCUGCAGAGAGGAAGCCCUUCUGCUUACAGGUUUUGGAAGGGCUUCCUCUCGAGGACCGCCUGAAU
Mature Sequence: CCUGCAGAGAGGAAGCCCUUC

Full insert: >MI0014186 TTGTTGTGCCCCGGGTAGCGCTCTGATGAGGGTGGAGCCTTCCCTTCCCAGGAACAGTCCTCGAGTTG
 GATTTGCCTTCTGCAGCCTATCTTCATGCGCAAGGACACCAAGATGAGTTCCAGTGGCGGATTCGAAAC
 CTCCCCATCCTAAGGATGTCTATAGTGTCTCTGTGGACCAGAAGGAGCGCTGCATCATTGTGCAGAACAA
 CCAACAAGAAGTGAGTAGCGTGAAGGCACCCTGTGTGCTTGTGTATGAGCGGGTGTAGGGTGGGGTTT
 GATCTGGCGTTAGGATCATTATTCAGGCCGGTCTGCAGAGAGGAAGCCCTTCTGCTTACAGGTATTGGA
 AGGGCTTCTCTGCAGGACCGGCCTGAATAATGTAATCAGCTCCATCGGGGAAGAGAACTTGAAATGA
 GTATCTAGTATAGTACCTGGCATAGATAGGGGCTCAGTAAATATTTGTTGAATGAATGAGATGATGAGGG
 AAATCACAGGACCATTAACCTTGGTCTCAGGTGGTTCTGCTGAAAGTTATTATGACCTTGGTCTCCTG
 ACTGGAGCCAGCCTTGGTAAGGGTGGCCTGGTGACATCACTGTCAGAGCTGGCAGGGAGCCATAGCTTG
 GAAGGCACCCCTGTCTGAAGTGGGAGTGACAGTCCTTTAAAGCAGTACTTG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)
Mammalian Cell Selection: Neomycin
Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq: [MI0014186](#)
Locus ID: 100422900



[View online »](#)