

Product datasheet for **SC401363**

MIR3163 Human MicroRNA Expression Plasmid (MI0014193)

Product data:

| | |
|---------------------------|---|
| Product Type: | MicroRNA Expression Plasmids |
| Species: | Human |
| Vector: | pCMV-MIR (PCMVIR) |
| Precursor: | UCCUCAUCUAUAAAAUGAGGGCAGUAAGACCUUCCUCCUUGUCUUACUACCCCAUUUUUAGAUGAGGAA |
| Mature Sequence: | UAUAAAAUGAGGGCAGUAAGAC |
| Full insert: | >MI0014193 ATGGTGAACCTCATCTCTACTAAAATACAAAAATTAGCTGGGTGTGGTGGGAGGCTGAGGTAGGAAAA TCACTTGAACCCAGGAGGTAGAGTTGCAGTGAGCCAAGATCGTGCCATTGCACCCAGCCTGGGCAACA GAGTGAGACTGTATCTCAAAAACAAAAACAAACACACATAGGGAATTACCGGGGAGTGGCATAGCAGAG TTGAGAGCATGGACCCTGGAGTCAGGCGGACAGGATGGAGTTCAGGCTTAGCCCTGATCAGCAAGTTCCT TAGCTTCTTCAAGCCTTGGTTTCCTCATCTATAAAATGAGGGCAGTAAGACCTTCCTTCCTTGTCTTACT ACCCCATTTTATAGATGAGGAAACCAACCAGTGGGGCTGTTGGGTAGAGATCAAATGCAATAACTTTT GTACAAATGAAATAACACACGTTTATTCTGTACAAATGAAATAACACATAAACAATGAAATCACACAGT CATGTGCCTGACACTTGGACACACCGTAATGGTGGCTCTGTTAATGGGACCAGTGATGAACAGTTAGAAG TGAACCCCTTGAGGGTAGAATATGTGAATCCGAGGCCCTGCAAGATGAGGGGACTCTCTCAGGTGGGTGG AGGAGCTGGTCCCTCGCATGTGTGATATTTGAGTTCGTAGGAGG |
| OTI Disclaimer: | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA. |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Mammalian Cell Selection: | Neomycin |
| Components: | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube |
| Also Contains: | Both mature and minor (star) miRNA can be over expressed from the construct. |
| RefSeq: | MI0014193 |
| Locus ID: | 100423029 |



[View online »](#)