

Product datasheet for **SC401377**

MIR3178 Human MicroRNA Expression Plasmid (MI0014212)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GAGGCUUGGGCGGGGCGCGGCCGGAUCGUCGAGAGCGUCCUGGCUGAUGACGGUCUCCCGUGCCACGCCCAACGCAGUCUC
Mature Sequence:	GGGGCGCGCCGGAUCG
Full insert:	>MI0014212 GGATGGGAGGGCTGGGGCGAGGCGGGGTCCGCAGAGATGTGGGGGGCGCAGAGATGGGGGGGTCCGCGG GCCCCGCGGGGACCGCGCAGCCTCGGGCTCTGTGCCCGGGGCCGCGGGTTAACTAGAGAAGGTGACGC GCGCCCTGTTGCTCAGCCGCGGCGCGTGCCTTCCCTGCGGCGATCCTAATTGGCCCGACAGAAGCCCGC CCCCAGAAGCCCTCCACAGTAGGCAAGGCCCGCGTCCGCAACGCTCCAGAGTTGCGCGTCAGCCT GGACTCAGCGCTGCGCTACGAGGCTGGGCGGGGCGCGCCGGATCGGTCGAGAGCGTCCTGGCTGATGA CGGTCTCCCGTGCCACGCCCCAAACGCAGTCTCTGGCTGGATGTGCTCGGAGCCACGAGGTGCAGAG GCTCGGACGGCGAAAGGGACCTGGCTGAGATCTGGGTCTTAACCGCGGGCTCGGCCAGCATCCCGGCA GCGTGCGCCCAAGTGACACCGCGCGGGCGCGGTCTGCGGGATCGCAGGCGCGGGTCTCATCTCCAG TTCACACGCGGACGAAACGGGCCACCCCTACCAGCTCTGAGAGCATGCACCTCCAACCTCCATCCCTCCC CACATGCCCGACCCCCAACCGTAACCTCAGGCAGACACCTGGGCGTCCCCCTT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0014212](#)

Locus ID: 100422974



[View online »](#)