

## Product datasheet for **SC401445**

### MIR3659 Human MicroRNA Expression Plasmid (MI0016060)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | MicroRNA Expression Plasmids   |
| Species:                  | Human  |
| Vector:                   | pCMV-MIR (PCMVMIR)   |
| Precursor:                | UCUACAAGCAGAUACAAGGAUGCCCUUGUACACAACACACGUGCUGCUUGUAUAGACAUGAGUGUUGUCUACGAGGGCAUCCUUGUGUCUGUGUGUGUG  |
| Mature Sequence:          | UGAGUGUUGUCUACGAGGGCA  |
| Full insert:              | >MI0016060 TCTGTGACTGTACATGGTAATTAATACATTAGCATTTATGCAAAATACCATCACAGCCGGGCTGGTATG<br>TCAGCTGCCTTCCTGAGCTGCAAAAATGTGAGCGCCCCGAGGTACCTCAGACAGGCCCCCGCCCTAGCA<br>GACAGCTCTGTTTGACCTGCATCATTCACTAACCTGTCTCGGGGAGCCTGTCTACATGCACACATCTGTG<br>TATGTTGGCCAATGTAATGCTGCTTCTTGTGAGTGTCTGTGGTATGCATATGTGCCTTGTGTGGAT<br>GCACACGTGCTGCTTGTCTCTACAAGCAGATACAAGGATGCCCTTGTACACAACACACGTGCTGCTTG<br>TATAGACATGAGTGTGTCTACGAGGGCATCCTTGTGTCTGTGTGTGCCTGTGTTGGTCTGGGCATGA<br>GTGCATATGAGTCTATGGTTTTATATGTGTACATCTCAGAATGTGTGTATGTACATTTTGTTTTTTTT<br>TCTTCCAGGTTGAAGCGATTCTCTTGCCCTCAGTCTCCCCAGTAGCTGGGATTACAAGCACCTGCCAACAC<br>GCCTGGCTAATTTTTAATTTTTAGTAGAGATGGGGTTTCACCATGTTGGCCAGTCTGGTCTCGAACTCC<br>TGACCTCAAATGATCCGCCACCTTGGCTCCCAAATGTTGGGATTACAGGCGTGAGTCACTGTACCT |
| OTI Disclaimer:           | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.  |
| E. coli Selection:        | Kanamycin (25 ug/mL)   |
| Mammalian Cell Selection: | Neomycin   |
| Components:               | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube  |
| Also Contains:            | Both mature and minor (star) miRNA can be over expressed from the construct.   |
| RefSeq:                   | <a href="#">MI0016060</a>  |
| Locus ID:                 | 100500801  |



[View online »](#)