

Product datasheet for **SC401452**

MIR3666 Human MicroRNA Expression Plasmid (MI0016067)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: AGUAAGGUCCGUCAGUUGUAAUAGAGACCCAGUGCAAGUGUAGAUGCCGACUCCGUGGCAGAGUUCAGCGUUUCACACUGCCUGGUCUCUGUCACUCUAUUGAAUU

Mature Sequence: CAGUGCAAGUGUAGAUGCCGA

Full insert: >MI0016067 ACTAATACTTTGTTTCATTTAAAGATACGATTTTTCTTCATATTTTATTACATCTTGTTCAAATTCGTG
TCATTGCTTTTTTCAGAGAGCTTACTAAAAATATTTTGTGTTTGACTTGTAATAATCAAGAACACAAAA
AGCACCATCACTGTTAAAGTTTAAAAAGCTTTTTATAAACCTGTTTACCATTTATAAAAGGTAGACTAAA
ACATAAGAGTTCAGCTGATGATACAGTCTTCGGTGTGACCAATTACGCCTATTTTTTTTTTTTCCTTCTG
CTGTGTCTGACTGATGAAAAGTAAGTCCGTCAGTTGTAATGAGACCCAGTCAAGTGTAGATGCCGAC
TCCGTGGCAGAGTTCAGCGTTTCACACTGCCTGGTCTCTGCTACTCTATTGAATTAGATTGATGATGGCA
GATTGCAGGGGGCACTAACTGGACCTCAGTGGGAATGCATGAAGTGTGTAGACAATCCTGGCAGGCACTT
TGCTTTGAGAACCCTTCACCCCTTCACAGCTGTTGGCACTAGTCCATTGCTAACAGTGTCCACAGGACTT
TAAAGAGACACCATGGGCCTTCTAATTTATGGTTTCAGTAACTTGTTACTGTGCGAGGCAACTGTGACCTC
AATCTCTTATTGACAAAAAGATGAAGTGTATTTGTTTTCTCTGACCTTTTAAACATAAGAAAAGCTATA
AAAGTTTTTCC

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms
should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0016067](#)

Locus ID: 100500896



[View online >](#)