

Product datasheet for **SC401481**

MIR3689d-2 Human MicroRNA Expression Plasmid (MI0016835)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	ACUGGGAGGUGUGAUCUCACACUCGCGUGGGAGGUGUGCUAUCGUCUCCUGGGAGGUGUGAUCUGUUCUCCUGAGCG
Mature Sequence:	GGGAGGUGUGAUCUCACACUCG
Full insert:	>MI0016835 GATCCCATGCCTCCTGGGAGGTGTGATATCATGGTTCCTGGGAGGTATGATATCGTGGTTCCTGGGAGGT GTGATCCCGTGCTCCCTGGGAGGTCTGATCCCGTGCTTCCTGGGAGGTATGATATCATGGTTCCTGGGAG GTATGATATCGTGGTTCCTGGGAGGTGTGATCCCGTGCTTCCTGGGAGGTGTGATCCTGTGCTTCCTGGG AGGTGTGATCCCATGCTTCCTGGGAGGTGTGATCCTGTGCTTCCTGGGAGGTGTGATACCATGCTTCCTG GGAGGTGTGATCCACGCTCACTGGGAGGTGTGATCTCACACTCGCTGGGAGGTGTGCTATCGTCTTCCC TGGGAGGTGTGATCCTGTTCTTCCTGAGCGGTGTGATCCTGTGCTCCCTGGGGGGTCTGATCCTGTGCTT CCTGGGAGGTGTGATATCATGGTTCCTGGGAGGTGTGATCCCGTGCTTCCTGGGAGGTGTGATATTGTGG TTCTGGGAGGTGTGATCCCGTGCTCCCTGGGAGGTGTGATCCCGTGCTCCCTGGGAGGTGTGATCCCGT GCTCCCTGGGAGGTGTGATCCCGTGCTCCCTGGGAGGTGTGATCCCGTGCTCCCTGTGAGGTCTGATCCT GTGCTTCCTGGGAGGTGTGATATCGTGGTTCCTGGGAGGTGTGATATCGT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0016835](#)

Locus ID: 100616344


[View online »](#)