

Product datasheet for SC401534

MIR3939 Human MicroRNA Expression Plasmid (MI0016596)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CUGGCUUCCAAAGGCCUCUGUGUGUCCUGUAUGUGGGCGUGCAGUACCUGUCACAUGUGUACGCGCAGACCACAGGAUGUCCACACUGGCUUCCAAACAUC
Mature Sequence:	UACGCGCAGACCACAGGAUGUC
Full insert:	>MI0016596 GGAGGGGAAGGTCGTGCGACCCGCTGGCCGGCGGAAAGGAGCAGCGGGAGCCGCGCCGTAGGGAAGCC ACGGGAAGACGGGCGGCTGCGCTGCGATGCGCGGCCCTCTCCCTGAAGCCCTCGGCTCAGGCGCCACAGG AAAGGCTCCTTCGGGCAGGCCCTCAGCTTTCAGCTTAAGTTTTTCTTTGGCCGAGTCAACTGCGGTGC CGAGTTGATTTTCTTTTCAGAGGTGTGTCGAGAACCATTCTTAGGAGACGCCTGCCGCGCTTTGGG ATAGGGCCTGTCTGCCATGTGGCTTCCAAAGGCCTCTGTGTTCCTGTATGTGGGCGTGACGTACC TGTCACATGTGTACGCGCAGACCACAGGATGTCCACACTGGCTTCCAAACACATCTCTGTGTTTCTGTCT GTGAGTGTGTACGCAGGTGTCAAATGTGTACGCACAGATCACAGGGTGTCCACACTGGCTTCCAAACACG TCTCTGGGTTCTTTATGTGCGTGTGTATGCATTGTGTGTACACACAGATCACAGGGCATCCACAC TGGCTTCCAAACATGTCTCTGTGTTTCTGTCTGTGGGTATGTACACACATATCATATGTGGAAGAGGTAG GAAGAGGCTGGAACAGGTAAGGGTCATGAAGGTGCAGTGCACATTGTACTGTTTGGCCTTTATATGTGT TTGGTT
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase .
Disclaimer:	The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0016596
Locus ID:	100500857



[View online »](#)