

Product datasheet for SC401551

MIR4253 Human MicroRNA Expression Plasmid (MI0015860)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVIR)
Precursor:	CCAGCCAUCGCCCCUAGAGGGGCCUAGGACUUACUUGUGCAGGGCAUGUCCAGGGGGUCCAGGUCUGC
Mature Sequence:	AGGGCAUGUCCAGGGGGU
Full insert:	>MI0015860 CACTGGGCAGCCTTTTCTGCCACATATACTACCTGTGTTGCTTCATCTCCTGCCAAGCCTCCCATCCCCT GTCCCTTTGTGAATCCAGACCTTGACTTTATCTTCCATGTGACTCAGCAGAAGGGGGTAAGGAGCCAGAT GCCAGCACTGTGGTGTCTGCAAGGTGGGAGCATGGCCACCCAGGATGCCCTTTCTTGAGCACTTTTCC ATTTCTCTTCCCACAGGAAAGACTGGCTGTAAACGGACGGACCCCAATCCTCACTCTGAGCTTTTCTG CAGGTAGACCTGGGACTCGGCCAGCCATCGCCCTTGAGGGGCCCTAGGACTTACTTGTGCAGGGCATGTC CAGGGGGTCCAGGTCTGCTCTGTAGTAGCCATTGCGGCAGACACAGTTGGTGGCCCTTCAGAAGTGGTC CGGCTGTTGATGGGACAGTGGGTACAGGCCTCATCCCCTTGGTTGGCCTTGAAAGTCCCAGATGGACAAC CTGAGAGAGAAGACAAGGAGAGTCACCAGGAACTGGATGAAGACTCAGAGAGGTCACGAGCCACGCAGC CTCACGGAGGGGCCACACAGCCTCGCGAGGGGCCATGCAGCACCCTGCTGAAAATCTCATCTCTTGA AGGGTTTGAAGTCTCCAGGCCCAACGTGGGAGTTCTA
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0015860
Locus ID:	100422914


[View online »](#)