

Product datasheet for SC401557

MIR4259 Human MicroRNA Expression Plasmid (MI0015858)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GAUGGGCCCCUUGUGUCCUGAAUUGGGUGGGGCGCUCUGAGUGGGGAAAGUGGGGGCCUAGGGGAGGUCACAGUUGGGUCUAGGGGUCAGGAGGGCCCAGGA
Mature Sequence:	CAGUUGGGUCUAGGGGUCAGGA
Full insert:	>MI0015858 CCAAAGCCAAAGTGGAGAGCAGGGCACGCCAGGGCGGGTCAATCCTTCAGAGCCGCGGTTGGGCTGGAGC GTGGAGCTGCAACGGTTGGAGCGACGGGCTCCGGTTCCTAAGCCACCAGCACAAAGGCAACTGCCCGCTGG AGGTGGGAGGGGCGGGGCCGGGGCGGGGCTGAGCCGAAGAACGCGGGCTAGTGGTTGCCAAGGTAACG CGTCAACACTAGGGCCTTGGCGGCGCCGAGGCGGCAGCAGAAGTCCGGAGTCAGGGCGTGTGGCTGAGGA GATGGTGAGGGAGAACCTCAGATGGGCCCCTTGTGTCTGAATTGGGTGGGGGCTCTGAGTGGGAAAAGT GGGGGCTAGGGGAGGTCACAGTTGGGTCTAGGGTCAGGAGGGCCCAGGAGTAAGGAGAACACCTCTAG AAGAAGCAAGGGAGATCCTGAGGTGGGTTCAAAGAGGAATTAGGGATAGCATCGATAACGATTGGAAT GTGGAGCGCGCTGAGCGGAGGAGTTGGGGAGAAAGGATGCCCAACCCTCTCCCCGTTGCCATTACAGC GGGGCCTCTGTGACCAAGATTGTGGATGTGGTGTGGGAAAAGCACTGGTGGGCCGGGCGCGGTGGCT TACGCTGTAATCCCAGCACTTTGGGAGGCAGAGGCGGGCGGATCATGAGGTCAGGAGATCGAGACCATC C
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in
Disclaimer:	miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0015858
Locus ID:	100422852


[View online »](#)