

Product datasheet for **SC401561**

MIR4263 Human MicroRNA Expression Plasmid (MI0015876)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: AUAGUGCUCUUCAGGGUUUUACUUGGGAGAUUGGAGUGGCCAGUGUCCUAAACAAUUCUAGUGCCUUGGCCACAACAUAC
Mature Sequence: AUUCUAGUGCCUUGGCC

Full insert: >MI0015876 ATATCCTTTGATTTTTCTGTGTGTGTTTCGTTTGCTTTTTCTTGTTATCCATCATGTGGTTCTGTTTT
 TGGTGTTCCTGGTAGTTTTGAGTGAATGTCAGACATTCCAGTTAAAAAACTACAGAGGCTCCAGATAACG
 CAGTCTTCTCCAGAGAGAGGGTTTACTGTCCCTCTGATAGGACATAGGCTGATGTGATTGGATGGTC
 ACTTTAATCCAGTCAGGGTTTTGCTGAGTCGAAGCTAGTTTATGGCTTTGGTAAGAGTCTGTTACTTC
 TGGTTCACCTTCTGTTCTGGGATAGTGCTCTTCAGGGTTTTACTTGGGAGATTGGAGTGGCCAGTGTTCCT
 AAACAATTCTAAGTGCCTTGGCCACAACATACTTCAGAATAGAATTCTCTGTCGTTTGTCAAATATTTA
 AACACACAAGCAGTGTGGATCTCAAATTCTAAATGTTGTCTTACAAGACTGCTGATTGCTTTGCTCTGCT
 TTTTCAGAGGCGTTTTCTTAATTTTTAGCTTCTGACCTGGGCTACGCCAGAATTCATCAAAGTCTTGAGAG
 GAGCTGGCTGCATGTTTGAGGTCTCTCAACACTCTGCCAAAAGCTCCAATGGTTCCTCTAAGCCTCAGCA
 GCAGTTTTGTGACTGTGCATAGCCTGGATTCTCAGCCTCCCTCCTGTACCCTA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0015876](#)

Locus ID: 100422965



[View online »](#)