

Product datasheet for **SC401608**

MIR4309 Human MicroRNA Expression Plasmid (MI0015837)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: UCUGGGGGUUCUGGAGUCUAGGAUCCAGGAUCUGGGUUUUGAGGUCUUGGGUUGUAGGGUCUGCGGUUUGAAGCCCCUCUUG
Mature Sequence: CUGGAGUCUAGGAUCCA

Full insert: >MI0015837 CATGGGGCTCCTGTGTCTCATCCATTGGTGGGATGCGGCTCCCTCCCGGGACCTCTCAGGGCACTGAT
 AAGAGCCACCAGGACTAGGTTTCAGCCACAGCAGGTGGGATGCACCCGCCCTTGCGGGCCGGGGTGG
 AGACCGGGAGTTTCTGGAGAAGGCCTTTCAGATTCTGAAGGCAGGCACTCCATGCCCCATCCCTGGAGCG
 CGAGTGACAAGTGTCTTTGTACAGAAGGACAGAAAGGCCCTACGTGGCCCTGAGCAGCTCTGGGGAAC
 ATTCTTGGTCCGGTTCCTGTCTGGGGTTCCTGGAGTCTAGGATTCCAGGATCTGGGTTTGGAGTCTTG
 GGTGTAGGGTCTGCGGTTTGAAGCCCTCTTGTAAACCAGGACAGAGGTGAAAGCAGTGGGGTGGAGTGG
 GGCAGCCCCCAAAGTTCAATGGACCTAAAGTGTGGTGTCTAGCCCTGCCAGGCAGCTGGCAGGTGGTC
 CAGAAACCCCTCGTTCTATGGAGACTGAGGCCCAAGTCTTTGGAGGTGGTCAGGGGTCTGTCCCAGAGCA
 CCTTTGGGCGAAAGGTAGGGGGAAGACCACCAAGACGACCCTCCCCAACACCAGCAAGGGAGCCAGCA
 TCTCCTCTGAGGCCACTGCCCTGTGGGTCCAGCCTCTCTGGGCTGGGTGTGTG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0015837](#)

Locus ID: 100422954


[View online »](#)