

Product datasheet for **SC401663**

MIR4444-1 Human MicroRNA Expression Plasmid (MI0016787)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVIR)
Precursor:	GUGACGACUGGCCCGCCUCUCCUCCUGGUCCCAUUAUUGAACUCGAGUUGGAAGAGGCGAGUCCGGUCUCAA
Mature Sequence:	CUCGAGUUGGAAGAGGCG
Full insert:	>MI0016787 AGAAGAAAGGGAAAGGCCTTTCTCTCTGGTGACGCGCATCACAAGGCGTTCCCTCTACGTGGACGCTCA CCATTGGCGGAGACCGCACTGACGTCCCGGAACTCGCCCTGACGCGAACGGCGGCAAAGGCCCTCTCCC CTGAAAGCTGGCGAGGCGCGCGCTCCCGCGGCACGCCGCGCTCCGCCCTCCGCCCGCAACCCGCTG GTCCAGTCCCTCTTCGCGCCCGTTGGCTCCGCCCCCTTCCCAATTGCTCCTCCCCAACAGACTGGCG GCGCGCGGAAACGCGTCACGTGACGACTGGCCCCGCTTCTCTCTCGGTCCCATATTGAACGAGTT GGAAGAGGCGAGTCCGGTCTCAAAATGGAGGTAAACCGCCGCCCGGTGCGCCCCAGCCCGACTCCGGCC GTCGCCGTGCGCGCCGGGGGAGGAGGTATTAGGGGAGAGCGGGGGTTGGTGGGAATGGCCGGCGTT GGGGGCTGGTTCGGTGGGAGCGGGGAGCGGGTGGACCGGTGCGCCGTCCCGCGCTCTTGCCTTGAG ACAGGCTGTGGAGGGGAGGGGAGCGAGAGTTGGAGACCTCCTTGGGGCCCTCGCCCGCTCGCCTTC ACCTTCGGCTGCACAGCGGCCCTTGGAGGCGCGCTTGGTCG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0016787
Synonyms:	MIR4444 hsa-mir-4444-1
Locus ID:	100616394


[View online »](#)