

Product datasheet for **SC401720**

MIR4498 Human MicroRNA Expression Plasmid (MI0016860)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMV MIR)
Precursor: AGGGCUGGGCUGGCAGGGCAAGUGCUGCAGAUUUUGUCUAAGCAGCCCCUGCCUUGGAUCUCCCA
Mature Sequence: UGGCUGGCAGGGCAAGUGCUG

Full insert: >MI0016860 ACCTAGAGAGAGGTGAGGAGGCAGGTCAGAGGCTCTTCTGAAAGATGATGCCTCACACGGTCCAATCCCA
CCCCTTCCAGGGGCTGATGGAATCCACATGGTATTGCCAGCACCTGATACTGATGAGAAGAATGGCCTG
AACCTTCTGCGGAGACTCTGGGTGGTCAAGTTTGACAAGGAGGAGGATCCGGAAGCTGGCTGAGAGGT
GAGAGAGCCAGCATGTCGCTTTCTCCTGCTGCATCCTCTCTTCTTTAGGAAATAACTGAACCCAGTGTG
CTCAGCCTCCAGTGGGAAGAAGGGCTGGGCTGGCAGGGCAAGTCTGCAGATCTTTGTCTAAGCAGCCCC
TGCCTTGGATCTCCCACAGGCTCTGGTCAATGATGGGCCTAGACCTGCAGCCAGACCTCTGCTCCTTGCT
GATTGACGACGTGATCTATCATGAGGCGGCTGTAAGGCAGGCAGGGCCGAAGCCCTCTCCAAGCAGTG
GCACGTTACCAGCGGCAGGCGGGGAGGTTATGGGCAGGCTCATGGAGATTTACCAGGAAAAGCTCTACG
TGAGTGACCTGTGCAGCCCTGATGCAGTCTGGGTGTGCGGTTGGGTGGTCTCTGCTACCCCCGCTAC
TGCTCCAGACAGGAATCTCGGGGCTCACGGTGGCTC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0016860](#)

Locus ID: 100616179



[View online »](#)