

Product datasheet for **SC401731**

MIR4508 Human MicroRNA Expression Plasmid (MI0016872)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVIR)
Precursor: AGGACCCAGCGGGCUGGGCGCGGAGCAGCGCUGGGUGCAGCGCCUGCGCCGGCAGCUGCAAGGGCCG
Mature Sequence: GCGGGCUGGGCGCGC
Full insert: >MI0016872 GGGTCCGAATGGAATGAGTGCAAGGCTGGGCGATGGCAGAGCTGAGGTGAAGTCCGGGCAGTCTAAGAAT
 ATCCTAAGAGGAGGGAGCAAGTTCAGAAAGTGAAGTGCAGGATGGAGAGAAGCTGCAAAGTCTGCAGTTT
 TTTACTCCAAACAGCTCAGGCATGCGCTGCCTAGCCCCGTGTGGTTAACATGTCCGTATTTAAAGAACAC
 ACAATTAATTGCCTAAAAGCATTATCCTTCTCAGCACTCTCCTATCCTTCAGATCTAACCTTTAGGCACT
 GCCGGCTGCCGGCGGGTCCAGGACCCAGCGGGCTGGGCGCGCGGAGCAGCGCTGGGTGCAGCGCCTGC
 GCCGGCAGCTGCAAGGGCCGAGAAATTAAACGCGAGTGTCCAGAATCGGGCATTGTGGGTTTAAATCCT
 GAATCCACTGGTCACTATCTGAGACTTTCCCCTGTACAAATTAGTTACCTTCTTTATCCCATCTGACA
 ACGGGAATAATAACAGCAAGAATTACACCGACTAGTTTACAGGTTTGTGTATTAAATTTATAATAA
 TTCATTAATTCTTATAAAGTGATTGGAACGATGCCTGATCTTTAAGATACCTGTAAGTGTCTCTTCTATT
 TTGAACCAACCAAGAATCCTGTGACTGAAGTCTTTGTG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0016872](#)

Locus ID: 100616275


[View online »](#)