

## Product datasheet for **SC401740**

### MIR4515 Human MicroRNA Expression Plasmid (MI0016881)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids  
**Species:** Human  
**Vector:** pCMV-MIR (PCMVMIR)  
**Precursor:** GCGGGAGGUGUAACAGGACUGGACUCCCGGCAGCCCCAGGGCAGGGGCGUGGGGAGCUGGUCCUAGCUCAGCGCUCCCGGA  
**Mature Sequence:** AGGACUGGACUCCCGGCAGCCC

**Full insert:** >MI0016881 GGGGAGCAGGGGCCCCAGAGAGACGGTGAGGGCGGGCGGGCGGGCCCGGGGGCCCCGGCTCCGCCTCA  
GCCCCCGACGCCTGCTCCCCAGCTGCGGCAGGCCCGAGTGAGGCCATCCTCCAGCTGCGCGGTTGCCAC  
GTTATGGACAAAACCTCCGCCGCATCGCCCAGGCCGCCTCCGGAGGCCGGCGCTGCCTCCCTGCCTTCCG  
GGAAAGGCGCTTCCGGGGGCCTCGCGCTCCGAGGCCCTCCCTGACCGCGCGCCCCGACGCCCCG  
GGTTGCCGCTGTGCGCGGCCGCGGGAGGTGTAACAGGACTGGACTCCCGGCAGCCCCAGGGCAGGGGCGT  
GGGGAGCTGGTCTAGCTCAGCGCTCCCGAGGTGGAGCTTCTGCAGGCTGGGTGCTGCCGGGCCTACGG  
GCACCCACCCCTGGCTCAGTGCCCTGGGTATTCTGACTCAGTAGGTCTGGGGTAGGGCCTGAGAATTTG  
TAATTTTAAGATGCATCGTGGATCTGCAGGCCGCTCTCAGGCTGCCTTTTTCTCTGCAGGCCCTTA  
TCACCAAGTCCAGCGCCCCATTGTGTATCTGGGACATGTCCAAGCCTGTCCCGAGGAGTGTTTTGCC  
TGGAAGTGGAGGGCAGGAAAGTATGGCCAAGGGAACACGTGAGGGAGGGG

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)

**Mammalian Cell Selection:** Neomycin

**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.

**RefSeq:** [MI0016881](#)

**Locus ID:** 100616404



[View online »](#)