

Product datasheet for **SC401764**

MIR4536-1 Human MicroRNA Expression Plasmid (MI0016906)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVIR)
Precursor:	AUGUGGUAGAUUAUUGCAGCAUUAUUAUACUGCCUGCUUUUAUACAUACAUAUACUACCUAUAUCGUGCAUUAUUCUACCACAU
Mature Sequence:	UGUGGUAGAUUAUUGCAGCAU
Full insert:	>MI0016906 GTGCGTCGGCTCCATTCCCACTGACGACGGGCTGTGTGTGTGTTGGGGTGGGAGGCCAGGGGTGGCGG GAGGGTGATTCTTGAAGGAGTCAGCTGTCTCCAGCCTTCCTTCTGTCTTGGCTTACACAGCGCGTT CATTCACTCGCCCTGTTTTAGAAGGCATTTGTGTGTATACTAATCAATTTGTGGTCCCATGTAGGACCC GGATAATATTGTTATGGCAGATGTCTGGATGATGTGAGGAGATGGGACTGCAGAAACAAGCAGGGTGAA AGAGGTGCAGAGCCTTTAAATGTGGTAGATATATGCACGATATATATACTGCCCTGCTTTTATACATAC ATACATACATACCTATATCGTGCATATATCTACCACATTATGTTTATACAGTGCATATATCTACCACA CATATATGTATGTATGTATGTATGAAAGCAGGGCAGTGACAGGATCAGAGGAATTGCATTTGTGCAGCCT CTCAGAAAAGCAAACTGTTTCCTGGCAGACCACATAGGAAGCAATTGGAATAAGAAGCACTGGGAACCTAA AGTATGCCCTTAAGGCTGTTGCCATGGGGATGGAGAGGCAGAAAGACATTTAAAGGTGACCGTTGAAGGG CTTTGCAGACAGAATGGATGAAGTGGGAATGTGGAGGGGATGATGGAGTTCAGATTG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0016906
Synonyms:	MIR4536 hsa-mir-4536-1
Locus ID:	100616155


[View online »](#)