

Product datasheet for SC401795

MIR4656 Human MicroRNA Expression Plasmid (MI0017284)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	AGGCUGGCGUGGGCUGAGGGCAGGAGGCCUGUGGCCGGUCCAGGCCUCCUGCUUCCUGGGCUCAGGCUCGGUUU
Mature Sequence:	UGGGCUGAGGGCAGGAGGCCUGU
Full insert:	>MI0017284 TGAGAAGAATGCCTGCAGCAGCGTGGGCACGGCCTGGGCACACTGGGCCACGCGGGCACAGCCGGCCATG GGCTGCAGCAGCTGGTGGAGTGGCGCCAACCTGAGGACACAGTGTGGCCGTGAGTGCTGCCCCGCTCCCCG ACCTGGTGACTGCAGGCCGGGGCCACGTCGCCACTGCTCCCTCTCGCACAGCTGAGCCCCGCTCCCATCCT CACTCATCGCCTCTCGGGACCTCCTGCTCCCCTGGGTTTGGGGAGACCCCTGGCGTCAGCACACCCGGGA CAGTAAAGCCATCGGGACAAAGGCTGGCGTGGGCTGAGGGCAGGAGGCCGTGTGGCCGGTCCCAGGCCTCC TGCTTCCTGGGCTCAGGCTCGGTTTGTCCATCTGCAAAATGGAGACCCACAGGGCCAACCTCCAAGCTT ATGAGGGTTGATGTGGGCAGCAGAGAGAGCTGAGACCAGGCCCGCCCCGGCAGACCCGTGGTGCAGAA GCAGGGCCTCTATAATGGTGATGACGGGGACGCAGCCTTGTGTGGGACACAAAGCACCCCTGGCCTCAGCC TAGGGCTGAGCAGTGAGCCCCACCCAGTGGCTGTGGCCGGGCTCTCTGAGCCACTGGTGGCCCCCTGGC CCTAGGGCCCCGTACCTCTCAGTGGCGCCACTGGCCTTGGGGCGC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0017284
Locus ID:	100616465


[View online »](#)