

Product datasheet for SC401799

MIR4659b Human MicroRNA Expression Plasmid (MI0017291)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVIR)
Precursor: CUGUUGACGUUGCCAUGUCUAAGAAGAAAAUUUUUCUCCAAAGUUUUCUUCUAGACAUGGCAGCUUCAGCAG
Mature Sequence: UUGCCAUGUCUAAGAAGAA
Full insert: >MI0017291 GGTTATAAAGGTATAAATAATGTCAAAAGTTATGTGAAAGTTAACAAATATGACAACTGACACAAATGAT
 ACGATGCTTTGGAGGTCAGTGGTTAAGTGAGAAAGTGGGATTCATTTCTGGAGAGTGGATGGAGTCAAACA
 CAAGAAAGCAAGGGGCCAGACTTCTCTTGCACAGCTATAGACTTAGTTATTTTCAAGATTGGAGTACAT
 TTGCAAGCTGGGAGAAACCCCTGCAGAGAAGTAAGGCAGATGACTGAGGCAACGATGGAGACAGGGTAGA
 CGGAATCAAGAGCTTAGAACTGTTGACGTTGCCATGTCTAAGAAGAAAATTTTCTCCAAAGTTTTCTT
 CTTAGACATGGCAGCTTCAGCAGTTTCTTCCAAAGAAATCCAAGTCAAGAGACTCCACATTAAGATTTTG
 AAATCAAAAGGGATATTGGGTGCACTCTTCCAGAGAAATATACTTTCTGAAAGCAAAAAATCCTTGAA
 TGTGAAGAGAACAAAAATCCGTGAAAACCATGAGGTAGTCAAAAGGAGAAAAATAAAGAAGTCTCCAC
 ACAGACTTGGCTCGGCTGAAATGAAATGGCCTGGCTTGGCTGCAGAAGGCAGTCAGCATTGATCTTACGA
 TGTCTCCACTGCTAAGGAAGGGAGAAAACAAACAGGGTTGCC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0017291](#)

Locus ID: 100616372


[View online »](#)