

## Product datasheet for **SC401805**

### MIR4663 Human MicroRNA Expression Plasmid (MI0017292)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | MicroRNA Expression Plasmids   |
| Species:                  | Human  |
| Vector:                   | pCMV-MIR (PCMV MIR)  |
| Precursor:                | CUGUGGUGGAGCUGAGCUCCAUGGACGUGCAGUGGCAUCUGUCAUUGCUGCCUCCUGGAGCUCAGGCCCUUGCAG  |
| Mature Sequence:          | AGCUGAGCUCCAUGGACGUGCAGU   |
| Full insert:              | >MI0017292 TTGACTCCTCAAGCTTGCCCTTGCTCAAGCCTTTGCCATCCAGTCACCTCCACTTTGAACATTCTCCCC<br>TGCACTCCCCAGCCCAGCCTTAAACGTCCCCTTCCTGTGGATCTCCACCCAAAGTCACTCTCAGACAA<br>GATGCTCCCGATCCCCAGACTGGCTCGGGTTCCCTTCAAATGTCATCATATGCCCTGTGCTTTTTTTTT<br>TTTAATGGCTTATTGGAATTGTAATTAATTTATTGTATAATTCCTCCTGGCCAGATCACAAGCTCTAGGG<br>CGGTGTCTCTGGTTTGCCTTCTGTGGTGGAGCTGAGCTCCATGGACGTGACGTGGCATCTGTCTATTGCTG<br>CCTTCCTGGAGCTCAGGCCCTTGCAAGGAGCTGAGGCATGGCAAGCCACAAGCGCACAGCCTGTCCC<br>CACCAAACCCCACTTCTGGTGCCAGAGGTGGCCCTGTACCCTGCACAGAACCATGTTTCAAAAAAATGTG<br>CTTCCTTTTAGGAAAACCTTGGGCAAAAAATCTCTACAAAACAATCTAAGGCCACCATTATTTAAAAGA<br>TGTTTTAAAGGCAGCTTATAGTGCTTTTGAACCTGAATCTAGGTAGTCTGTTTACACCCTGCTTTCCA<br>GGGATTAATCAAAAACAACAAAACAATGAGCAAACAGAGGGATG |
| OTI Disclaimer:           | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.  |
| E. coli Selection:        | Kanamycin (25 ug/mL)   |
| Mammalian Cell Selection: | Neomycin   |
| Components:               | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube  |
| Also Contains:            | Both mature and minor (star) miRNA can be over expressed from the construct.   |
| RefSeq:                   | <a href="#">MI0017292</a>  |
| Locus ID:                 | 100616260  |



[View online »](#)