

## Product datasheet for **SC401807**

### MIR4665 Human MicroRNA Expression Plasmid (MI0017295)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids  
**Species:** Human  
**Vector:** pCMV-MIR (PCMVIR)  
**Precursor:** CUCGAGGUGCUGGGGACGCGUGAGCGGAGCCGCUUCCACGGCUCGCGCCGCGGCGGUAGCCCCGCCACAUCGGG  
**Mature Sequence:** CUGGGGACGCGUGAGCGCGAGC

**Full insert:** >MI0017295 CGCCAGGATGCGGTAGCCCTGCTGTAGCTCGTAGGTCAGCTCCTGCTCCTTGCAGCAACCGCCTCCGAT  
CACCATCGCCTCCATCTCTTCTCCTGATCGTCCGCGTCTCCAGCGAGGAGGCACTCCTTCCGTGGGCC  
GGCCCTGAGGTCTGGGCCGCGCTGCCACCTCCTCCTCGTCTCCTCCTTTCGGCCGCGGTGGCGGCC  
GCTCTTCTCCCCAGCCGGCTCCATCGCTCCCGGCGTCCCGGGCACACTCATGCCCGGCGAGGCCTAGGC  
TGGGCGGTGTGGAACAGCCGCTCGAGGTGCTGGGGACGCGTGAGCGCGAGCCGCTTCTCACGGCTCGG  
CCGCGGCGCGTAGCCCCGCCACATCGGGCCTGGCTCTCTCAGCCGCGGCTCGGCTCGGTGCGCGCGG  
GGTCTTGAGCTACCGCCGGCGGGGCGGGTTACATGGCGCGGAGGAGGGGGAGAGAAGAGGAGAGC  
CTAGGTGCCCTCTCTCCCTCCCCCGGCGGCGCGCGCACGGCCCTCGCCGGCCTACCCCTCCG  
CACCCCGCCGCTGCTCTTACCGCCAGGCCCTGCCGGTTCGCCGCGGCTCGGCGCCCCACCCCTC  
CCGGCCAGCGGGGGCTGGGAGGGGGCGAGGAGGAAGGGGTTCTCT

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)

**Mammalian Cell Selection:** Neomycin

**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.

**RefSeq:** [MI0017295](#)

**Locus ID:** 100616288



[View online »](#)