

Product datasheet for **SC401808**

MIR4666a Human MicroRNA Expression Plasmid (MI0017296)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	AUCACUUAAAUACAUGUCAGAUUGUAUGCCUACAAAUCCCUCCAGACUGGCAUACAAUCUGACAUGUAUUUAAGAGAU
Mature Sequence:	AUACAUGUCAGAUUGUAUGCC
Full insert:	>MI0017296 AGGCAAAGAAGGTAAAACAGGATCCTTGCTACTAAGATTTTATAATGGAGAAAATTCAGAGCAATCGTGT CCAGAACCCATTTCTCTAATTATAGTGTTACAACAGAGAAAGACAAATGTTAATGATAAATGAAGGCAGT GCACATCAAATATGTGAGAATCCACAGAGACTAAGGGAATTAGGACTCACTCAACAAAATGGAAGCATT TCTGGGTTGGAGAGACATGCCAAGATTCTTAACTAAGCCAAATTTTAACTAAGCTTCATCTCACCTA TCTGAGATGGGTGAGTGACAATCACTTAAATACATGTCAGATTGTATGCCTACAAAATCCCTCCAGACTG GCATACAATCTGACATGTATTTAAGAGATTGCCCTCACCCATCTCAGATGGAATTTCAAGTCAGATT GCTAATTATTATGCAGTAAATAAAATAATTTGGAAGGAATCATCTCCAATTAAACCACTTTGATTAAAA AACAAAAATAACAAATGTTGGTGAGAACGCAGGAAAGGGGATGCTTGTACACTGCTGGTAAGAATGTA AATCAGTACAGCCACTATAGAGACAAATATGGGCATTCTTTAAAAAATAAAAAACAGACTTACCATATGA TCCAGCAACCCATTGCTGGGTATATATCCAAAAGAAAGGAAATCAGTA
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0017296
Synonyms:	MIR4666 hsa-mir-4666a
Locus ID:	100616308


[View online »](#)