

Product datasheet for SC401815

MIR4672 Human MicroRNA Expression Plasmid (MI0017303)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: GGCUGCUUCUCGCCUCUGUCCAGCUGUGUGGCCUUGGACAAGCCUCUUGGUUACACAGCUGGACAGAGGCACGAAACAGCC
Mature Sequence: UUACACAGCUGGACAGAGGCA
Full insert: >MI0017303 TCAGCAGCTGATGGATGATGATGGGGGTGGGTGACCTTGAAGTTGAGGTCAGACTGCTGCATCCAGGCC

CAGCCTTTGCCCCATGGTAATCTGTAAGGTAGAAGCAGCCGCTGGGTGGTACGGAGTCATCTCCCTGAA
 GGCAGCTGACGCTGTACCCCTGAATCCTTAACACGAGGAGCTAGGGATGAGCGCTAGGATGGGTATCG
 CAGGCGAGGACAAAGGGCTTGCCCCAGGTACACAGCTGAACAGAGGCAGGGGTGAGCCTGGGTCCAGGT
 CAGTCTGACGCCGGAGGTGAGCTGCTTCTCGCCTCTGTCCAGCTGTGTGGCCTTGGACAAGCCTCTTGG
 TTACACAGCTGGACAGAGGCACGAAACAGCCTGACCTGAGAGGCAGTGCCCCAAGGCCTTGACAGGGGTC
 CAGGGCCAAAGCTGCGCAGAAACCTTGGGGGCTGGCCAGGCAGTTGGAAGCAGGAGACAGGGTTTGAA
 GGTCCCAATCCCATTCCTAAAGCTGTGGCAGTCATGGGTGACTCTGGGTGGCTCTGCCCTCCCAACCAG
 GCAGCTGTGGTAGCTGGATTGAGCGAGGCAGGGGCTGCAGAACCTTATCCTCTTGCCCTTCCAAGGCGT
 CCCATGGATGACATGGATGACATGCTGTTCCAATGACACGGACATTGCTCA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0017303](#)

Locus ID: 100616429


[View online »](#)