

## Product datasheet for **SC401818**

### MIR4675 Human MicroRNA Expression Plasmid (MI0017306)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids  
**Species:** Human  
**Vector:** pCMV-MIR (PCMVIR)  
**Precursor:** CAUGAGAAAUCCUGCUGUCAACCAUAGCCCUGGUCAGACUCUCCGGGGCUGUGAUUGACCAGCAGGACUUCUCAUG  
**Mature Sequence:** GGGGCUGUGAUUGACCAGCAGG  
**Full insert:**  
>MI0017306 GTAAGAGAACTATTGGTTCCAAGAGACAGAGACAAGGGATTTTGTACACATTTTGGGAAAAAATCAA  
AATGAGAAAGATACCTGCCCCCTCCCATGGCAACCCTGCCACCTGATATTATGTGCAAAAAAGCAGCA  
TGGCATCTAATCTAGTGTTCATAGGGCTTCTAGGGACAATAGCTCCACCTCACAGGAGTTCAGATTCA  
ACTGAAAGGCATGCAGGAGTTAGATTTATGAAGATTAAGTCAAGCTTTATAACTGATAAAATCTAGGC  
TGAATTGGACAGACTCACTCCATGAGAAATCCTGCTGGTCAACCATAGCCCTGGTCAGACTCTCCGGGGC  
TGTGATTGACCAGCAGGACTTCTCATGCACTGAGTCTTAGAGCCCTCTCCAAGCAAAGGAGAGTAGAAG  
CGACTTACTCAACCATGGCATAGTTCCTCCCATGTCAGGAAAACCACTTACCTGGAGAGAAATAGGGTGT  
GGGCAAGGAGCTTTACCTTCTCCAGGAGAAATCCTGCAAAGGAAAGGAAATCCCGGCTACTCACCTT  
GGAATCAAGTGGGGAACACATTCTCCATTTTCAGTGTTCTTTCTGGGAAATGCACAATACCTTTCAAGG  
CCACAATCCACATCTGTTTGAGACTGGAGCCCTTATTTTTGGCATT

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)  
**Mammalian Cell Selection:** Neomycin  
**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube  
**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.  
**RefSeq:** [MI0017306](#)  
**Locus ID:** 100616383



[View online »](#)