

Product datasheet for **SC401863**

MIR4720 Human MicroRNA Expression Plasmid (MI0017355)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVIR)
Precursor:	AAGCCUGGCAUUAUUUGGUUAUAACUUAAGCACCAGGUAAAAUCUGGUGCUUAAGUUGUACCAAGUAUAGCCAAGUUU
Mature Sequence:	CCUGGCAUUAUUUGGUUAUAACUU
Full insert:	>MI0017355 TACTGAAAATCATGAATGATGTGATTTGTTCCATGTATCTGTGGAACATCCCTCACCTTCTCCGGTATG CGCCTTCTGCACTCCACACGTGGTCGCCCTCCATCCTCCCTAGTGTCTGTGAGCTTGCCTGGTCACCAAGT CAGGAGGAAGAAGATGAGGAGGCTCACCAGTTGTAGGCAAGACTCTTCTCCTCGGGAGATCTTGTAAATA TGGCCAGTTTCAGGGCACAGTTGTACACATTCTACTTAAATGGAGGAGTTCAGTGATGCCGTGGCCAAA AGAGTTATGCCATTTCTAAAAGCCTGGCATATTTGGTATAACTTAAGCACCAGGTAAAACTGGTGCTT AAGTTGTACCAAGTATAGCCAAGTTAACTGTCGTTGTGATTAATAAGATTTAATTTTTATGCTTGTGC TTTTCAGAGTAAAGGGCTTGCTTTCAATAGTAGAATGACAAAGCTTCTCCTCCAATGACAATCTGTTAGC ACCCTGTTGACTTAGCTCAGGACATGAGACTGACTTAGCTCAGGACATGAGACATTTCTGGAAGTGATT TCCATGAGAACAAGGCTTATTCCTTTAACTTAAACCCAAACTTTGTTATTTTAATTATTTTCAGGCAA TCATCTTTATATACTTTTCCTTACAAACCCAGTGTTGAGAGCATAG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0017355
Locus ID:	100616150



[View online »](#)