

Product datasheet for **SC401894**

MIR4751 Human MicroRNA Expression Plasmid (MI0017390)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CCCGGAGCCAGAGGACCCGUAGCUGCUAGAAGGGCAGGGGUGUGGCUUCUGGGGCUGGUCUUCAGCUCUGGCG
Mature Sequence:	AGAGGACCCGUAGCUGCUAGAAGG
Full insert:	>MI0017390 AAGTGGGGATGCCGCCTCTGCCCCGCCACAGCAGCCCCCTCCTCTTCTCCACCTCAACCTTCTCGCCT GGCCCCCTACCCACATCCTGCCACCACCCGAGGGGACCGCAAGCAAAAGAAGAGAGACCAGAACAAGTCG GCGGCTCTGAGGTACCGCCAGCGGAAGCGGGCAGAGGGTGAGGCCCTGGAGGGCGAGTGCCAGGGGCTGG AGGCACGGAATCGCGAGCTGAAGGAACGGGCAGAGTCCGTGGAGCGCGAGATCCAGTACGTCAAGGACCT GCTCATCGAGGTTTACAAGGCCCGGAGCCAGAGGACCCGTAGCTGCTAGAAGGGCAGGGGTGTGGCTTCT GGGGGCTGGTCTTCAGCTCTGGCGCCTTCATCCCCCTGCCTCTACCTTCATTCCAACCCCTCTCGGCCG GGTGCAGTGGCTTATGCTTGTATCCCAGCACTTTGGGAGGCCAAGGCAGGAGGATCGTTTGAGGCCAGG AGGTCAATACCAGCCTGGGCAACATAGTAAGACCCTGTCTCTATTAATAAAAAAAAAAATCAACCTTCTTC CCCACCAACCACCAACTCCTCTCTACTCTTATCCTTTATCCTCTGTCTCTGCTTATCACCTCTCTTG CGTATTCTGGATCTCCTTCCCTCCTTCTCGTCCAAATCATGA
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0017390
Locus ID:	100616483


[View online »](#)