

## Product datasheet for **SC401910**

### MIR4767 Human MicroRNA Expression Plasmid (MI0017408)

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | MicroRNA Expression Plasmids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Species:                  | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Vector:                   | pCMV-MIR (PCMVIR)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Precursor:                | ACAUGGCCCCGCGGGCGCUCCUGGCCGCCGCCGACUUCGGGGCCAGCCGGGGGAGAGCGCGGGAGCCCCGAGCGU                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Mature Sequence:          | CGCGGGCGCUCCUGGCCGCCGCC                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Full insert:              | >MI0017408 TACCATCAAAGTGCGCCACAGCAGCAACCTTGGAGCAGGCAGGCAAGTATCCTCCAGGCTACTTTGCAG<br>ACGGGGAAGCTGAGGCTGGGCACCGCTCAAGGTCGCACGCATGCGTAACCGAAAAGAGAGCGCGGAGGTC<br>CCGCCGCGCTCGCTCTGAGGAGAGGACGCGGAGCCGCGCGCCCGCCCGTGGACCCAGGGACTGCTCC<br>CACGTGCTCCGGCCCCACGCGCGGGCGCCGACCGGAAGCGGGTGTGGGGCCGCGTAGACGCCGCGGCC<br>ACGCGCGCCCCGAGCCCGGACATGGGCCCGCGGGCGCTCCTGGCCGCCCGCCGACTTCGGGGCCAGCCG<br>GGGGCAGAGCGCGCGGGAGCCCGAGCGTCCCTGCATGAACACCGCCCCGCGCGGCCCGCCAGGCCGTGAC<br>GTACCCGCGCGGACCGTCCCCACGCCACACAAGACCGCCCTTACTGGAGGCGGCGGCTGCACACTACC<br>CACCCAGAAGAAGTCCGTCCATGTCAAAGATGAGGTGGGTGACGGGCTGCGGGGGCGCCGCCATGGTGGC<br>GCCTTCTGGGTCTGGGTGGGGGCGAGGAGGAAGTGCAGCGCACCCGCCCGCCCGCGCACGCGCACTGA<br>GGACCTTCTCGCGCGCCCCGGGTCTGCTGCGCACGCGCTTCCGCGG |
| OTI Disclaimer:           | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Components:               | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Also Contains:            | Both mature and minor (star) miRNA can be over expressed from the construct.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| RefSeq:                   | <a href="#">MI0017408</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Locus ID:                 | 100616467                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |



[View online »](#)