

Product datasheet for **SC401922**

MIR4776-2 Human MicroRNA Expression Plasmid (MI0017420)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: CUUAUAGCAGUGGACCAGGAUGGCAAGGGCUCUACUGUCCUUUCAGGAGAGCCCUUGCCAUCUGGUCCACUGCAUUAUG
Mature Sequence: GUGGACCAGGAUGGCAAGGGCU

Full insert: >MI0017420 TAAATTCTCACTATTTCATTGGTTTGAATGCTTGTGTATATCCCAAGGGTATACACAATGACCTCTACACA
 CATGACCTCTGTTTTACATAGTCAACTTTTATTTGCTGAGGCAGAAATAATGCATGCAGCAAAGACTGC
 CAGCTGTTTACCAAAATCCCTTCTCCTCTTCTCCTGAGCACACAGCTAGACTATACTTCCCAGTCTCCA
 TTGCAATTCTGTGTGGCCACATACTAGAGCTCTTGACTGCTGTGTTTTGCATCCCCTTCTACCTTCTG
 GCTTCCTGAGATGGAGCAGACTATATGCAGTGGACCAGGATGGCAAGGGCTCTACTGTCCTTTCAGGAGA
 GCCCTTGCCATCCTGGTCCACTGCATATAGCCGAGCTCCCACTGACCTGGAACCTCATCCTGGACTGTT
 GTCCAAAAGAGAGATAGATGTGTTTTCCCTCAAACATGAGAAGCTGAATTTGGCGGTCTATTTGTTACG
 ATATCATAGCCTGTCTAACTGAAACAGGGCACATGTAGTATGTAGCACAAAGTTTCTTAGCTAGTGAGAG
 TGCTTAGCAGAAGCTCAGTGCCTGCCCTTACCAAGATTTGAGGTTTCCTACTCATGGAGTGTCATCAGT
 AACTCTTGCCAAGTGATAAATCTTTTGAATATTGGTGAATATGATAACA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0017420](#)

Locus ID: 100616472



[View online »](#)