

Product datasheet for SC401938

MIR4792 Human MicroRNA Expression Plasmid (MI0017439)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GCAGCCCGUGAGCGCUCGCGUGCCUGGCAGUGCGUGGAAGAACAGGGCGGGUGGGGCCGCGCACAUUCUCG
Mature Sequence:	CGGUGAGCGCUCGCGUGGC
Full insert:	>MI0017439 GCTCCCGGCGACCGGCGGCCGCGACACGGTGGACCCGGGGCTCGGTGAAGCCCTGGCTGTGCCCCGGG TGGCCGAGTAAAGGGCGCCACGGTGAGGAGACGACGGCCGGCGCCACGGCTTGAGGGGCTGCCTTAG GCTGCAGGACTGAGACTCGAGCTGCTCGCTCGGGCTTTACACCTTCGGTGTCTTCTGCTGAAGCTCCC TGAGCCCGGCCCTGGCTCTACGAGCCGGCTGAGAGCGCTGGGCGGCCCTTCCGGGTCTCGGGACTCGAG TCCCCTCCGCTGCGCGTTAGGCAGCCCGGTGAGCGCTCGCTGGCTGGCAGTGCCTCGGAAGAACAGGGC GGGTGGGGCCGCGCACATCTCTGCAAGAGTGAGCTCCTAACTCTGGGGCTCCAAGCTCCCAGCAAGTCTT ATTTAATGGAAATGTTACTGAGACTGGCTACACTGTCTGGGGAAAACCCCATGCAATTGTTGTAATGTA TATTTTCATTGCAGTATCAAAGGTGCCTGGTGCACAGCTAGATGCTTGTTCACGAATTCACACACCTGTT TAACTAGCAGTCACACCAAGCAACAGAACATGACCTGCACTCTCAATGCCCCCTCAGCTCTTCCACCAA TCGGGGCAATCCCTACCCTGCCTTCCAACACCAGGACTAGTTTT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0017439
Synonyms:	miR-4792
Locus ID:	100616448


[View online »](#)