

Product datasheet for SC401946

MIR4800 Human MicroRNA Expression Plasmid (MI0017448)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: GGAGAAAGGAGUGGACCGAGGAAGGAAGGAAGGCAAGGCUGUCGUCCAUCCGUCCGUCGUCCACCUACCUGUCAGUCC
Mature Sequence: AGUGGACCGAGGAAGGAAGGA
Full insert: >MI0017448 ACAGGCCCACTGCTGTGCCATTCTGGAAGCTCCAGCTGCTGCTGGGCTGCCTGGCACTGCCCCGCTTGCCG

GTCAGGGCCTGCCGAGCTGCCTGCCCCCTCCAGCTGGGCAGAGTCCCCTGCAAGGAGGCAGGGCCAGCT
 TCCACATCCGGAGCCCTGGTCAGCATAGCCGCCACGGTCTGTTCTCAGATTCCATAATCATTCCAGAAGT
 ATTAACGTCATTGCTGCAAACCTCGGCAGGTGCCGTGTGAGGGGCTTAATGACCACCACAGGGAGCTCA
 GACCCCAACCTGGATCCCAGGAGAAAGGAGTGGACCGAGGAAGGAAGGAAGGCAAGGCTGTCTGTCCAT
 CCGTCCGTCTGTCCACCTACCTGTGAGTCCACATAGGCTCCTGGCGTGGACAAGGGGTCTGTGAAGGGCG
 GGAAGTGGGTGAGCACCTGGGGCAGGTGGGTGGTTAAGGTCCTTCCCACTTCGCAGGTGTCAGAACCTAG
 GGCTGGGCTCTCGGGCCAGGCAGGCCAGCCAGCCACGCCAGCTGGGCAGCGTCTGCTCTGGTAGGA
 CTGTCAGACGCACACGCGCACGACCTAGACACACCCACTCATGTACATGCTCACACATGCAGACACACC
 TGGCGTCCCAGGTCACATGTTCTGGGGATGATGGCCTTCAGGGGTCAT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0017448](#)

Locus ID: 100616358


[View online »](#)