

Product datasheet for SC401956

MIR5002 Human MicroRNA Expression Plasmid (MI0017868)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVIR)

Precursor: UCUUCCUCUCUGUCCUCUGGAAUUUGGUUUCUGAGGCACUUAGUAGGUGAUAGCAUGACUGACUGCCUCACUGACCACUCCAGAUGAGGGUACUC

Mature Sequence: AAUUUGGUUUCUGAGGCACUUAGU

Full insert: >MI0017868 GGAGAAGGTGACTTTGTTGTGTCTGGGAAGGAACCAGGGCCTTGTAACTCAGAGAATTGGGGAAAAGCTT
 TGTGGCACCTTCATTTTTCTGGCTCCTTCATGGCTGTTTTGGCAGGGATTTTTGCCCCGTCTCACATTCC
 AAGTTCTTTCATGAGAGCTTTGTGGTCAGAGGCAGGGAAGTACAGACCTTAACCTCTGAAGACACA
 AATTCTCCCTGGCAGCCCTGTCTCTTTGGGAGTTGTTCTGGTTCCTCTGCTGGGGCATTTACACTGT
 TTTCTCTGAAGTCAAATCCATCTTCTCTGTCTCTGGAATTTGGTTTCTGAGGCACTTAGTAGGTGA
 TAGCATGACTGACTGCCTCACTGACCACTTCCAGATGAGGGTTACTCTCTCACAGCACGGCCTTGAGGC
 AGCTCATGCCAGGGATGCTGCCTGTGTTCCAAGTGGATCTGCATTCTGGGAACAGCCAATTGATCTGA
 TGGTGCGGGGTTACTTTTTATTTTTTAAATTTAAAAATTAATAATATCCTTAGTGGTGATTAATA
 CTTCAGAACTGTTTTGAATACTAATTTGGTTTTTGGAAACAGCCCACTCCATTCGGGAAAATCTGGGG
 AAGAGGGAATCAAATCTGGCAAGCTGGGTGACAGTATTTAGTTTAAAAACAAAACAACCTT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0017868](#)



[View online »](#)