

Product datasheet for **SC401970**

MIR5090 Human MicroRNA Expression Plasmid (MI0017979)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: UCUGAGGUACCCGGGCGAGAUUGGUGUAGGGUGCAAAGCCUGCCCGCCCCUAAGCCUUCUGCCCCAACUCCAGCCUGUCAGGA
Mature Sequence: CCGGGGCAGAUUGGUGUAGGGUG

Full insert: >MI0017979 CATCCTAGAGGGGGCTAGGATGCCCCGAGTCCTAAAGTAGTTGCAGCTTTTTTTTTTTTTTTTTTTTTT
 TTTTTTTTTTTTTTTGTTAAGTGGTGTGGAGTCGGGAGTAGTGGCGCGCGCTTCTAGTCCCAGCTACTCGA
 GAGGCTGAGGCGGGAGGATCCCTCAGCCCAGGAGTTGGAGGCAGCTGTGAGCTATGATCGCGCCTAGGCG
 ACTGAGTGACACCCCGTCTCTAAAAATAAGTAACGGGGGGATGGGCGGGGCGGCTACACTTGTACGAGA
 AATGTCAGGTGAAAAAGCCTTCTGAGGTACCCGGGCGAGATTGGTGTAGGGTGCAAAGCCTGCCCGCCCC
 CTAAGCCTTCTGCCCCAACTCCAGCCTGTCAGGATTGGAGCTGCTTCCGAGCACCTGGACCCCAAACCT
 CCTGTGCCCGCTGACGCAGCTGCAGGAGCTTGACCTGTCTAACAACCACCTGGAGACGCTGCCGGAAC
 CTGGGCCTGTCCACCTGCGTGTCTCCGCTGCGCAACAACAGCTGGGGGATGTTACTGCCTTGTGCC
 AGTTCCCAAGCTCGAGGAACCTCAGCCTGGAGGGCAACCCCTTCTGACGGTAAGTGGGAGCCTCCACC
 AGCTTCATACCTGGGGCCAGGACTCTGTTGGCTCAGGCTCAGCATCTCCTGAGC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0017979](#)



[View online »](#)