

Product datasheet for **SC401979**

MIR5186 Human MicroRNA Expression Plasmid (MI0018165)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: UCAGCCAGCUUUGACUUGACCCUCACCCUGAUUUUCUACCAACCUUUCUCAGCUGAUUUUCUUGGGGAGAGAUUGGUAGAAAUCAGGUGAGAGGGUCAUGCC

Mature Sequence: AGAGAUUGGUAGAAAUCAGGU

Full insert: >MI0018165 TTCTACTATTCATCTCAAATCTTTCCTAATCTATCATTGTGGGGTAGGGTTTCACTTCCCTGCTTCTA
ATTCATAAACTATCATGTTTATCATTATTTGGTCAGATTTCTGATTTTCTAGCTGATACCTCCCATTCTA
ATCCAAATTTGCTTTCTCCCTCTGCCAGGGAGTTCTCATGACATGTAGCTAGGGCCTCACTTTGTGGCA
ATTTTATAAGCCTCTCATTAAACATTGGCACTGCCTCAGTTGCTTTCTGCAGTATTCTCCAGGACTGG
ACCCCTGCCTACTCCCAACTCAGCCAGCTTATGACTTGACCCTCTCACCTGATTTCTACCAACCTTTCC
TCAGCTGATTTCTTTCTGGGGAGAGATTGGTAGAAATCAGGTGAGAGGGTCATGCCATAAGCTGGCTAAC
ACTTTCATGCATTAATTAATACAGGAATCATAGCTATCATCTTAATCTATCACTCTCGTATCTATCTGAT
ACTATGCATTAACGTTCTCAAAGGCAAATCCAACCTAGAACAGTGTTTCTCAAATACCAATCAGTGGGC
CAGAAGCAACAGAATTAGCTCAGTAACTACTAGAACTCCAGGGTAGGTGCTTTAGTAGGTCTAAAAATA
TAAAGCAGAACATTTCTATGTTTGTGTTTTATAACTTCCCCCAGGTGAGTCTGATGCACAAAAAAGT
TTTGAACCCCTGGAACAGAT

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect
the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0018165](#)



[View online »](#)