

## Product datasheet for **SC402012**

### MIR548ar Human MicroRNA Expression Plasmid (MI0019131)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids  
**Species:** Human  
**Vector:** pCMV-MIR (PCMV MIR)  
**Precursor:** AAAAGUAAUUGCAGUUUUUGCUGUUGAACGUAGUGGUAACUGCAGUUUUUUUGC  
**Mature Sequence:** AAAAGUAAUUGCAGUUUUUGC

**Full insert:** >MI0019131 GATATACATGTTGCTCTTACTGATCTTGATAATTTGCAAGTTTAAAAATCCATTTTGGCACATTGTTTT  
GTAAAAATTTATTCTCAATTATGAGACTCACACTGTAGATACTAGGTGTTAATATGCTGCACGCAGTAGC  
TAAGCCATCTCCACATGTATCAGTTTTCCAGTGTACCTAATTTAAACTTTAGCTAAATGACATATGCTA  
TTTTTTGTAGATCTTGGTTCCGTTTAAAATCCTTGAAATTACAAAAATTTAATTGGCAAATACACATAA  
CCTATTATTAGTTTGGTGCAAAGTAATTGCAGTTTTTGCTGTTGAACGTAGTGGTAAACTGCAGTTA  
TTTTTGCACCAACCTAATAGGATTTCTGTTGGTTTTGATTCTTTTAAATAAGCGGCCTACTTATTTTAT  
GTATTACAATTCTGTTGACTTAAAGTAAATGTATGTGATAATGTTTATTTTTATTCCTTTCAAATAAACC  
AGAATTTTTCTTTCAATTTGTAACCTAATTATGGCCTTAACTGAATCAACCTGAAGATAATGGAGTGAAG  
AAAGTGAAGTGTACAACCCACAGCATCATATGTGATACCATAATGACTGTCTACACATAGCCGATCGTC  
GTGAGTGCTGTACCTGCTGGTTTGGG

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)

**Mammalian Cell Selection:** Neomycin

**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.

**RefSeq:** [MI0019131](#)



[View online »](#)