

## Product datasheet for SC402022

### MIR548h-5 Human MicroRNA Expression Plasmid (MI0016751)

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | MicroRNA Expression Plasmids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Species:                  | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Vector:                   | pCMV-MIR (PCMVMIR)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Precursor:                | ACAAAAGUAAUCGCGGUUUUUGUCAUUUUUAAACUGUAAAAACCACGGUUGCUUUUGC                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Mature Sequence:          | AAAAGUAAUCGCGGUUUUUGUC                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Full insert:              | <pre>&gt;MI0016751 AGTACTTTTTGTATGTAATAAAAGAAACCAACTAAAATTAGATCATGACAAAAAGGAGGGACATTGTTTT AAGAATAAAGAAATATTTAAACCCAAAGACAGTAGGTGAGGGAGAGTAGAGCCACGCAGGTGGGGGCTA GGGCCCACAATCAGGTCACTTCATGTTCCCTCCTGTCTCCGTTTCTTTCCAAAAGCTCGCTCCGTGCTT GCCTTTCTGTGTATGGCAGATTCCTCTGCTGACTCACCTTGCGCACAGCACAGCATGGTCTCCCTAAGCC CCTTTTCTATTAGGTTGGCAGAAAAGTAATCGCGTTTTTGTCACTTTTAACTGTAAAAACCAAGG TTGCTTTTGGCCAGCCATAACGTATGATCTTCCAGCTATGCTGCTGCAGAAATCTCAAGCAACAGA AAACCTCAAACCAAGCAAGGAAACCTCAGTCCCAGCTGAGTCAATGGGTTGCTGGATTCTGTGAGGCC TCCCTCTCTGTTCAAATTAGTCATGGGGCAGTGCCATGTCCCAGGATCTGTGGGACTTCCCCTTCTAGTG GCTGTGAGGAGAATAGGTTTTCCAAAGTGGACACAGGATAGGTGTAACAATGATGGTCAAAAGGACTAT GGCTCTTAATTTGGAATTGGTACCTTTCCG</pre> |
| OTI Disclaimer:           | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.                                                                                                                                                                                                                                                                                                                                                                                                  |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Components:               | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Also Contains:            | Both mature and minor (star) miRNA can be over expressed from the construct.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| RefSeq:                   | <a href="#">MI0016751</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Locus ID:                 | 100616455                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |


[View online »](#)