

Product datasheet for **SC402033**

MIR548z Human MicroRNA Expression Plasmid (MI0016688)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVIR)

Precursor: AAGUUAUUAAGUUGGUGCAAAAGUAAUUGAGAUUUUUGCUACUGAAAGUAAUGGCAAAACCGCAAUUACUUUUGCACCAACCUAAUAGAUGCCAAUG

Mature Sequence: CAAAACCGCAAUUACUUUUGCA

Full insert: >MI0016688 GAGGTTGCAGTGAGCCAAGATCGTGCCACTGTACTACAGCCTGGGCAACAGAGCAAGACTCCATCTCAAA
 AAAAAAAAAAGAAAGAAAAAGAAAATTTTCACAGTAAGCTTGTGTTGAATTAAGAAATTGGGGATTAAGGAG
 GGCTCATAATCTATCAATCATAAAGTCAGAAGTAACTCTTCACATCTTTTAAAAAGATTCACTTGATAAA
 AAGATTCTGAAAGAAAAATTTAAAGTAGGAGTTAATCAGAATAAAAGAGATTAAGCAACAATTTTCTAC
 TTAGGCTTCAACAAGCCTGAAAGTATTAAGTTGGTGCAAAAGTAATTGAGATTTTGTCTACTGAAAGTAA
 TGGCAAAACCGCAATTACTTTTGCACCAACCTAATAGATGCCAATGCCGATTCTTTCTCATCGCAATT
 TTTCAACAAAAGACAGTAAGCCTATAAGTAAACAGTATAAATTATCTTTGGTGGGTAAGGTTTCCAGCAG
 TTTGTCCTACCATCTTTTTCATACCTAGTATACTTACAAAATAAGTTAGGACTGGCCAGGCACGGTG
 GCTCACACCTGTAATCCCAACACTTTGGGAGGCCAAGGCAGGCGGATAACGAGCTCAGGAGTTTGACACC
 AGCCTGACCAACATGGTGAACCTGACTCTACTAAAAATACAAAAAAAAAAAAAAAAATAGCTGGA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0016688](#)

Locus ID: 100500856



[View online »](#)