

## Product datasheet for **SC402045**

### MIR5584 Human MicroRNA Expression Plasmid (MI0019141)

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | MicroRNA Expression Plasmids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Species:                  | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Vector:                   | pCMV-MIR (PCMVIR)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Precursor:                | CAGGGAAUGGGAAGAACUAGAUUUGAAUCCAGACCUUUGUUCUCCCUUGCCCAU                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Mature Sequence:          | CAGGGAAUGGGAAGAACUAGA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Full insert:              | >MI0019141 AGTCACAAAGGGAACCTCGGAGATGTGCAGGACAGAGAGAAGAGCTGACGGGACATCTGGTACCTTCGGTT<br>GGAGGGAGGGCAGTTTCTGAGGATTAAGAAGCTTTATTGAAGGGAGGTTTGAAGGAGAAGAATTGTG<br>GAGAAGAAAGCTACTGTGATTTTCATGCTGATTCCATGCCATTAACCTTCAGGCACTTTACATCCATTAT<br>CTAATTTTATCCTCAGGCACTCAGCAAGGCAGTTCACGTTATTTCCACCTTGAGATTCAGAGGAGGAAGT<br>GACTAATCCAAGATCCATAACAGGGAATGGGAAGAACTAGATTTGAATCCAGACCTTTAGTTCTTCCCT<br>TTGCCCAATTATGCTCCTTGTTGGTTTGTCTGGAGAAAGCCTGGCAGAGGGTGCAAGGGCAACAGTGGGGT<br>GGGGTGCAGAGGCTAGGAGACCTGACCAGAGTGTAGAGGGGTTACAGGGCCAAGGAGAGGAGATATCT<br>AGGGGGCTAAAGATGGGTGAGACCTCTGAGCCTCCCTAAACAGCATACTGCAACCTCTTCATCAGAAA<br>CCCAATTTTTTTTTTTTTTTTGGAGATGGAGTCTCACTCTTGTGCCCAGGCTGGAGTGCAATGGCGCAA<br>CTCGGCTCACTGCAACCTCCACCTCCTGGG |
| OTI Disclaimer:           | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.                                                                                                                                                                                                                                                                                                                                                                                                                      |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Components:               | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Also Contains:            | Both mature and minor (star) miRNA can be over expressed from the construct.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| RefSeq:                   | <a href="#">MI0019141</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |


[View online »](#)