

## Product datasheet for **SC402053**

### MIR5680 Human MicroRNA Expression Plasmid (MI0019280)

#### Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GCAUUGGGUUAGCAGGUUAGCCAGCAUUUCCUCCUGGACACACAGGAGGAGAAAUGCUGGACUAAUCUGCUAAUCCAAUGC
Mature Sequence:	GAGAAAUGCUGGACUAAUCUGC
Full insert:	>MI0019280 TCTCAAGAGCAACTCAAATTCATCAATTTAAAAAGTTGCTGCTTAAAAATAACAAGGCTGCAAAGGCACC GTGACCTGGCCTTCTCGTTGTAATGTTTGCAAAATCTTTAACATGACTGGTTATAGAAAGGACAGTTA CATTGAGGTCCAGACAGCTGCCAACTCTGATTTAAACACCAAGTAAAAGTTTCATCTAATCAAAAAGCACC CTTTCCTAGCTCTACCGGTTGCCTTCAATTGTA AAAACAAAAATATCTAGGAGGCTGTTAAACCATCATC AAATACTAAAACCTTTGCTTGCATTGGGTTAGCAGGTTAGCCAGCATTTCCCTTCTGGACACACAGGA GGAGAAATGCTGGACTAATCTGCTAATCCAATGCAAGCAAAAATTTGCATTGGAAAGTTGGATTTTTTGG TTTTTGGTTTTTTCTTTTTCTAATACTAGGCCCGACCCGAGAGAATCTAATCTAATTGTTAAAAAAA AAAAAATCCGGGTTCCAGGGTGTCCAGCTTTGTATGGGTATTTTTTTCACATGGGCCCCAGGAGACTC TGATGTGTGGCCAGGTTGAAACCACCATAGGATTAGCCAATCATCAAGCAGGCTAAGGGAACGTCC CGTTGATACAGCAAAAGGTCTGCACAGGTGTTGCAAACTCAAATGTCCCTGAAG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	<a href="#">MI0019280</a>



[View online »](#)