

Product datasheet for **SC402084**

MIR5706 Human MicroRNA Expression Plasmid (MI0019314)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: AGCUAGGUCUUCUGGAUAAACUGCUGAAGCUUCUACGUCAUUCAGCACUUGCUUCAGCAUGUUUCCAGAGGAUCUAGCU

Mature Sequence: UUCUGGAUAAACUGCUGAAGCU

Full insert: >MI0019314 TAAATGCAAAGGACAAGGTCCAAAGGAACTAAACGTTCTCCATTGAACACATAAATGATAAGAAAGCAA
GGTAGCCTTATTGCTTATGTGGAGAAAGTTTGAGTGGTCTAGATGGATCGAACTAGCCACAAGATTCTCT
TATGCCAAAGCCTAATCCAGAAACAAGGCTGTAACCTCGTCTTAATTTTCATGAAAGCTGAGAGAAGTGAGGA
AGCTGCACAATTAAAGTTTGAAGCTAGCAGAGGTTGGTCCATGAGGTTTAAGGAAAGGAGCTTTCTCCGT
AAACTACAAGTACAAGGTGAAGCTAGGTCTTCTGGATAACATGCTGAAGCTTCTACGTCAATTCAGCACTT
GCTTCAGCATGTTTTCCAGAGGATCTAGCTAAGATAATGGATGAAGGTGGCTACACAAAAGGACAGATTT
TTAAAGTACATGGAACAGCCTTCTGTTGGAAGAAGATGCCATCTAGGACTTTTCGTAGTTAGAGAGCAGTC
AGTGCCTGGCTTCAAAGAACAGACTGACTCTTGTTAGGGGCTAATGCAGCTGGTAAGTGAAGCCA
ATGCCCATTTACCAATTCTGAAAATCCTAGAGTGTTTAAGAATTATGTGAAACCTATTCTGCCTGCGCTCT
ACACATGGAATAGCAAAGCCCAATGGCAGCACCTCTGTTTATAGCATAAT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0019314](#)



[View online »](#)