

Product datasheet for **SC402120**

MIR6130 Human MicroRNA Expression Plasmid (MI0021275)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: CAAUGAAUGUGGAAUUGCAGGCAUCCCUUCUAUAUACUGAUGUAAUUUUUCCUUGGAUUUAUAUGCAGUACUGAGGGAGUGGAUUGUAUGGUAGUUCUAUUUUUAC.

Mature UGAGGGAGUGGAUUGUAUG

Sequence:

Full insert: >MI0021275 CCCAATCCCTCTCCTCTACTTAGCTCAGCCTCTGGTAACCACTACTCAACTCTCTGCTTCTACGAGTTC
CACTTTTATGATTCCACATATGTGAGATCATGTGGTATTTCTGCTGGCTTATTTCTCTTAGCATGAAA
TCATCCAGGTTTATTGATGTTGTACAAATGGTGGGATTTTGTCTGTTTGTGACCGAATACTCTTCCA
TTATGTAAATATACCACATTTTCTTTAGTGATTCTTCTGCTAATAGACACTTAGGTTGATTTTCATGCTT
GGTTATTGTGAATAATGCTGCAATGAATGTGGAATGCAGGCATCCCTTCTATATACTGATGTAATTTTC
CTTGGATTTATATGCAGTACTGAGGGAGTGGATTGTATGGTAGTTCTATTTTACATTTTCTGAGGAAAC
TCCATGCTGTTTTCCCTAATAGCACTAATTTATATCCCACTAACAGTGACAATTTTGGCACATTATTTA
CAATATTTATCTTCTGCCTTTTGTAAATAGTAATCCTAGTAGAAATTAGGTAATATTTTATTTTGGTTT
TGATTTTCCATATCCTAATGATTATGATGTTGAGCATATTTTAAATACAGGTTGGCCATTTGTATGTTT
TCATTTGAAACGTATCTATTCAAGTCTTTAGCCATTTTAAAAGGGTTATTTGTTTCTTGCTATTAAG
TTCCTTATA

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms.
This should not affect the function of the mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0021275](#)



[View online »](#)