

Product datasheet for **SC402158**

MIR6720 Human MicroRNA Expression Plasmid (MI0022555)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Human

Vector: pCMV-MIR (PCMVMIR)

Precursor: UUGAGCGAGAGAUUGUGGCGACCGAGUUCUCCAGCCUGGUAGGCGCCGCGGAAGAAGGGGAAGCGCGCCUGCAGGAACUGGUAGAUCUCGCUGAG

Mature Sequence: UUCCAGCCUGGUAGGCGCCGCG

Full insert: >MI0022555 GCGTCGTAGCCCGCGGGCATCATGTCGAGGCCGCGTAGCCGCCCTGGCTGTGGCAGCCGAGCGGCGCCG
ACGGGGGCGCCTGGAAGTCGAAGCCCTGGGGCAGCAGCGACGCCCCGAAGCCCAAGCCGCTCACCACGCG
GTGGTACATGGGCTTGAGCGCCTGGCACTTCCGCCTGAAGCCGCGCGGCCGCGGCGGAACGAGCCCTCC
TCGAACATGAAGTCTGCTGGCGGGTTCGATGGTCCAGTAGTGCCCTTGCCGGGCGCCCGAGGCCCTTAG
GCAGCTTGATGAAGCACTCGTTGAGCGAGAGATTGTGGCGACCGAGTTCTTCCAGCCCTGGTAGGCGCC
GCGGAAGAAGGGGAAGCGCGCCTGCAGGAAGTGGTAGATCTCGCTGAGCGTCAGGCGCTTGCTGGGCGAG
CTCTGGATGGCCATGACGATGAGCGCGATGTACGAGTAGGGCGGCTTCTCGGGCGCCCGCAGCCCCGAGC
TCGCTTCTTGGCGCCCCGCTCCCGCGCCCCGCGCCGCCCGCCCGCGCTTTCAGGCGAGCCGAGGG
GGCGCTGGCCGAATTGGAGGAGGACGAGGACGAGGCGCAGGAGGCGGAGGACGACGAGGAGGAGGAG
GTGGTCTCCGGGGCGCGCGCGCGCGCGCGCGCGCGCGCGCGCGCGCTCATCAGGGCGGCTGGAG

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in
Disclaimer: [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0022555](#)



[View online »](#)